3BU1-4U-JUN06-FDA-3

## **Premier Wrist Blood Pressure Monitor**

with Microlife Averaging Mode

Instruction Manual  
Model #BP3BU1-4U  
ID #3BU1-4U

PLEASE NOTE:  
THIS MEDICAL INSTRUMENT MUST BE USED  
ACCORDING TO INSTRUCTIONS TO  
ENSURE ACCURATE READINGS.

**microlife**<sup>®</sup>

# Premier Wrist Blood Pressure Monitor

Instruction Manual

Table of contents

---

## **1. Introduction**

---

- 1.1. Features of the BP3BU1-4U Blood Pressure Monitor
- 1.2. Important information about self-measurement

## **2. Important information on the subject of blood pressure and its measurement**

---

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?
- 2.4. MAM Technology

## **3. The various components of the blood pressure monitor**

---

## **4. Putting the blood pressure monitor into operation**

---

- 4.1. Inserting the batteries
- 4.2. Sleep mode
- 4.3. Setting the time and date
- 4.4. User selection

## **5. Carrying out a measurement**

---

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Discontinuing a measurement
- 5.6. Memory – storage and recall of the measurements
- 5.7. Memory – recall of 3 intermittent average results
- 5.8. Memory – cancellation of all measurements

## **6. PC-Link functions**

---

6.1. Installation and Data Transmission

6.2. Change user 1/user 2 memory during PC connection

## **7. Error messages/malfunctions**

---

## **8. Care and maintenance**

---

## **9. Warranty**

---

## **10. Reference to standards**

---

## **11. Technical specifications**

---

## **12. How to contact us**

---

## 1. Introduction

---

### 1.1. Features of the BP3BU1-4U Blood Pressure Monitor

This Microlife blood pressure monitor (with PC Link Connection) is a fully automatic, digital blood pressure measuring device for use on the wrist. This monitor enables very fast and reliable measurement of your systolic and diastolic blood pressure as well as pulse by way of the oscillometric method of measuring. This device offers a very high and clinically tested measurement accuracy and is designed for maximum user-friendliness. Before using it, please read through this instruction manual carefully and keep it in a safe place. For further questions on the subject of blood pressure and its measurement, please contact your doctor.

#### **Attention!**

### 1.2. Important information about self-measurement

Self-measurement means control, not diagnosis or treatment. Your values must always be discussed with your physician. **You should never alter the dosages of any medication without direction from your physician.** The pulse display is not suitable for checking the functioning of heart pacemakers!!

**In cases of irregular heartbeat, measurements made with this instrument should only be evaluated after consultation with your doctor.**

**Electromagnetic Interference:** This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens, etc.). These can lead to temporary impairment of the monitor's measuring accuracy.

## 2. Important information about blood pressure and its measurement

---

### 2.1 How does high/low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood

vessels is altered. Blood vessel width is affected by way of fine muscles in the blood vessel walls. Your level of arterial blood pressure changes periodically during heart activity. During the "blood ejection" (Systole), the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value). Your blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

## 2.2 Which values are normal?

Blood pressure is too high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, while at rest. In this case, please consult your physician immediately. Long-term values at this level endanger your health through continual damage to the blood vessels in you body. Should your systolic blood pressure values lie between 140 mmHg and 160mmHg and/or your diastolic blood pressure values lie between 90 mmHg and 95mmHg. Consult your physician. Regular self-checks will be necessary. If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician. Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately. If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician. Which values are normal?**

The following standards for assessing high blood pressure (without reagent to age) have been established by the U.S. National Institutes of Health JNC7, 2003.

Category	Systolic (mmHg)	Diastolic (mmHg)
<b>Normal</b> Self-Check	<120	<80
<b>Pre-Hypertension</b>	120 - 139	80 - 89
<b>Stage 1</b> Hypertension	140 - 159	90 - 99
<b>Stage 2</b> Hypertension	>160	>100

### **Further information**

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called "labile hypertension." Consult your physician.
- Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.

## **2.3 What can be done if regular high/low values are obtained?**

- a) Consult your physician.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constrictions caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. Secondary hypertension can be ascribed to specific organ malfunctions. Please consult your physician for information about the possible origins of your own increased blood pressure values.
- d) There are measures you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

### **1) Eating habits:**

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

## **2) Previous illnesses**

- Consistently follow all medical instructions for treating illness such as: diabetes (diabetes mellitus or sugar diabetes), fat metabolism disorder, and gout.

## **3) Habits**

- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea and chocolate).

## **4) Your physical condition**


- After a medical examination, and with your doctor's approval and direction; exercise.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your physician before beginning your exercise routine. You must receive advice regarding the type and extent of exercise that is appropriate for you.

## **2.4 MAM Technology**

Microlife Average Mode (MAM) Technology is a new concept for optimum reliability in self-measurement of blood pressure. A single measurement is rarely your real result. An average of three readings is medically proven to provide the most accurate result. With the push of a single button, the Microlife Average Mode blood pressure monitor takes three consecutive measurements, and calculates the average.

### 3. Components of your blood pressure monitor


The illustration shows your blood pressure monitor, consisting of:

### 4. Putting your blood pressure monitor into operation

#### 4.1. Inserting the batteries

After you have unpacked your device, insert the batteries. The battery compartment is located on the left lateral side of the device (see illustration).

- Remove cover as illustrated.
- Insert the batteries (2 x AAA 1.5V), observing the indicated polarity.
- If the battery warning appears in the display, the batteries are exhausted and must be replaced by new ones.

## Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AAA» Long-Life or Alkaline 1.5V Batteries. The use of rechargeable batteries is not recommended.
- If the blood pressure monitor is left unused for long periods, please remove the batteries from the device.


## Functional check:

Press and hold the ON/OFF button to test all the display elements. When functioning correctly all segments appear.

## 4.2. Sleep mode

The device is designed to stay at the sleep mode after batteries are inserted. In this mode, the time and user icon will be displayed. You can conduct the following operations in this mode.


- Reading the set date: Press the TIME button. The date set will be displayed for 2 seconds.
- Recalling your measurements: Press the MEMORY button. The measuring data stored can be displayed. For more information refer to the section on Memory storage and recall.
- Changing the User: (See 4.4 for more details) This advanced monitor tracks memory for two users separately. Press and hold the time button (⌚) until the 1 or 2 icon flashes. Next release the button. Press the memory button (M). The unit will now record readings for User 1 or User 2 as indicated.


### 4.3. Setting the time and date

This blood pressure monitor incorporates an integrated clock with date display. The advantage is the exact time each measurement is recorded. After new batteries have been inserted, the clock begins to run from the following setting: 2003-01-01 00:00 O'clock. You must then re-enter the date and current time. For this, please proceed as follows (Example: Entering 2004-06-20 Time 09:30 o'clock):


- 1) Press and hold the Time button for three seconds. User 1 or 2 will begin flashing. Press the Time button again briefly. Now you can begin changing the time and date.


- 2) The correct year can be entered by pressing the MEMORY button. (Example: 1 x press)

- 3) Press the TIME button again. The display now switches to the current date, during which the first character (month) blinks.


- 4) The corresponding month can now be entered by pressing the MEMORY button. (Example: 6 x presses)

5) Press the TIME button again. The last two characters (day) are now blinking


6) The corresponding day can now be entered by pressing the MEMORY button. (Example: 20 x presses)

7) Press the TIME button again. The display now switches to the current time, during which the first character (Hour) blinks


8) The corresponding hour can now be entered by pressing the MEMORY button. (Example: 9 x presses)


9) Press the TIME button again. The last two characters (Minutes) now blink

10) The exact time can now be entered by pressing the MEMORY button. (Example: 30 x presses)

11) After all settings have been made, press the TIME button once again. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.




## Further Information

With each press of the button (TIME, MEMORY) one input is made (e.g. switching over from hours to minutes mode, or altering the value by +1). However, if you keep the respective button depressed, you can switch more quickly to the desired mode, or find the desired value respectively.

### 4.4. User selection:

This advanced blood pressure monitor allows you to track blood pressure readings for 2 individuals independently.

- Before measurement make sure you set the unit for the intended user. The unit can track results for 2 individuals. (User 1, User 2)
- The unit is set to User 1. Push and hold the TIME button for 4 seconds to initiate the user selection.
- When the user icon is flashing, press the MEMORY button to toggle between users.
- When the correct user is selected wait three seconds for the user icon to stop flashing. The monitor is now ready for use.
- Repeat the procedure again to change users again.

## 5. Carrying out a measurement

---

### 5.1. Before the measurement:

- Avoid eating, smoking as well as all forms of exertion directly before your measurement. These factors influence the measurement result. Relax by sitting in an armchair in a quiet atmosphere for about ten minutes before your measurement.
- Always measure on the same wrist (normally left).
- Compare readings at the same time of day, since your blood pressure changes during the course of the day.

### 5.2. Common sources of error:

#### Note:

Comparable blood pressure measurements always require the same conditions! These are normally always quiet conditions.

- All efforts by the patient to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the wrist artery lies considerably lower (higher) than the heart, an erroneously higher (lower) blood pressure will be measured! (Each 15cm/6" difference in height results in a measurement error of 10mmHg!)
- A loose cuff causes false measurement values.
- With repeated measurements, blood can accumulate in the respective arm, which can lead to false results. Correctly executed blood pressure measurements should therefore first be repeated after a 1 minute pause or after the arm has been held up in order to allow the accumulated blood to flow away.  
(Does not apply to a MAM measurement.)


### 5.3. Fitting the cuff

- a) Remove all objects and jewelry (e.g. wristwatch) from the wrist in use. Put the cuff over the wrist with the monitor face

on the inside of the wrist.

- b) The distance between the cuff and the hand should be approx. 10 mm (1/2").
- c) Secure the cuff with the Velcro fastener, so that it lies comfortably and not too tight. No space should remain between the cuff and the wrist.

- d) Lay your arm on a table, with the palm upwards.

Support the arm a little with a cushion, so that the cuff rests at about the same height as the heart. Remain still for 2 minutes, before beginning the measurement.


#### 5.4. Measuring procedure

After the cuff has been appropriately positioned, and you've selected User 1 or User 2, the measurement can begin:

- a) If you want to take an Average Mode measurement position the slide switch to the 3. Remember, in this mode, the unit takes 3 measurements (one after another) and calculates the result. There is a 15 second rest time between measurements. The unit beeps five seconds before each upcoming measurement. (Data obtained in Average Mode are indicated by the MAM symbol.)

- b) If you want to take a Single measurement, position the slide switch to the left. In this mode there is only one measurement cycle.
- c) Once you've selected the correct user, and measurement mode you can press the ON/OFF button to begin. After pressing the ON/OFF button the pump begins to inflate the cuff. In the display, the increasing cuff-pressure is continually displayed.
- d) After reaching the inflation pressure, the pump stops and the pressure slowly falls. The cuff pressure (large characters) is displayed during the measurement. When the device has detected the pulse, the heart symbol in the display begins to blink and a beep tone is audible for every pulse beat.
- e) When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values as well as the frequency are now displayed. (Pulse=number of heartbeats each minute) Example (Fig.): Systolic 118, Diastolic 73, Pulse 75. The measurement results are displayed until you switch the device off. If no button is pressed for 5 minutes, the device switches off automatically to save the batteries.

f) When the unit is set to the Average Mode setting 3 separate measurements will take place in succession followed by a calculation of your detected blood pressure value. There will be 15 seconds of resting time in-between each measurement. A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings begin. If one of the measurements causes an error message, it may be repeated one more time. In case any additional error occurs, the measurement will be discontinued and an error code displayed.

## 5.5. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the ON/OFF button can be pressed at any time. The device then immediately lowers the cuff pressure automatically and enters into sleep mode.

## 5.6. Memory – storage and recall of the measurements

This blood pressure monitor automatically stores the last 30 measurements values for each user. By pressing the MEMORY button, average of existing memory readings displayed first with an “A” and (MR) at the lower left. Press MEMORY button again, the last measurement (MR30) as well as the previous 29 measurements (MR29, MR28, .....MR1) can be displayed one after the other. Data obtained in Average Mode are indicated by a MAM symbol.

Further

(MR1: Oldest measurement.)  
For User 1

(MR2: Values of measurement after MR1.)

## Information

When the unit is in Average Mode, time between measurements is carefully controlled to ensure accuracy. If you take successive measurements when the unit is in the single measurement mode wait for several minutes in a relaxed position before you repeat a measurement. Make sure all measurements are taken in a seated position.

## 5.7. Memory-recall of 3 intermittent average results

The blood pressure monitor can display the 3 intermittent average results of the last stored memory. By pressing the MEMORY button at least 3 seconds, the LCD will show  icon first with one short

beep. Then release the MEMORY button, and the 3 intermittent average results will appear. NOTE: This only occurs if the last stored measurement is taken in Average Mode.

## 5.8.Memory– cancellation of all measurements

### Attention!

Before proceeding to clear all the memory, make sure that you don't have to make reference to the readings in the future.

To delete all stored readings, press the MEMORY button for at least 7 seconds. The display will show the symbol "CL" and 3 short beeps will sound to indicate deletion of stored readings.

### Note:

If your last stored measurement was taken in the Average Mode, you will see  icon with a short beep during the deletion procedure. (The display will show these symbols during deletion of all stored measurements if the last measurement is taken from the Average Mode)

## 6. PC-LINK FUNCTIONS (PC-Link Software can be purchased as an optional accessory by calling 1-800-568-4147 or visiting [www.microlifeusa.com](http://www.microlifeusa.com))



This unit can be used in connection with your personal computer (PC) running the Microlife Blood Pressure Analyzer (BPA) software. Your PC will allow a capacity of monitoring 80 patients, each with 1000 data (note: overuse will lower system efficiency). The memory data can be transferred to the PC by connecting the monitor via cable with the PC

### System Requirements for Blood Pressure Analyzer Software:

- Windows 98 SE, 2000, XP
- CD-Rom Drive
- Minimum 256MB RAM
- 500MB Available Hard Disk Space
- USB Port version 1.0 or higher
- Microlife Blood Pressure Monitor with USB Port

### 6.1. Installation and Data Transmission

- Insert CD into CD ROM drive of your PC. The installation will start automatically. If not, click on SETUP.EXE
- Connect the monitor via cable with the PC. Three horizontal bars will appear on the display and last for 3 seconds.
- The bars will then flash to indicate that the connection between computer and device is successfully made. As long as the cable is plugged, the bars will keep flashing and the


buttons are disabled.

During the connection, the device is completely controlled by the computer. Please refer to the 'help' file for software instructions.


## 6.2. Change user 1/user 2 memory during PC connection

While your monitor is connected with PC you can switch between user 1/user 2 memory data by pressing the TIME key. Please note that new patient files need to be prepared on the PC for each user by entering names and other information.

## 7. ERROR MESSAGES/TROUBLESHOOTING

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).

Error No.	Possible cause(s)
ERR 1	No pulse was detected. Ensure the cuff is properly placed on the inside of your wrist and repeat a measurement
ERR 2	Unnatural pressure impulses influenced the measurement result. Reason: The arm was moved during the Measurement (Artefact). Repeat, keeping still and quiet.
ERR 3	Inflation of the cuff takes too long. The cuff is not correctly seated. Re-position cuff and repeat
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.

## Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains blank when the instrument is switched on.	<ol style="list-style-type: none"><li>1. Check batteries for correct polarity.</li><li>2. If the display is unusual, re-insert batteries or exchange them.</li></ol>
The device frequently fails to measure blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none"><li>1. Check the positioning of the cuff.</li><li>2. Measure blood pressure again in peace and quiet carefully following the details in Section 5.</li></ol>
Every measurement produces a different value.	<ol style="list-style-type: none"><li>1. Please read the following information listed in Section 5.2 "<b>Common sources of error.</b>" Repeat the measurement. <b>Please note: Blood pressure fluctuates continually so successive measurements will show some variability.</b></li></ol>
Blood pressure values measured at your doctors office differ from your measurements.	<ol style="list-style-type: none"><li>1. Record the daily development of the measured values and consult your doctor. <b>Please note: Individuals visiting their doctor frequently experience anxiety which can result in a higher reading than at home under resting conditions.</b></li></ol>

## 8. Care and maintenance, recalibration

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.

c) Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soap-suds. The cuff must not be washed!

d) Do not drop the instrument or treat it roughly in anyway. Avoid strong vibrations.

e) Never open the device! Otherwise the manufacturers warranty becomes invalid!

### Further information

Blood pressure is subject to fluctuations even in healthy people.

**Comparable measurements always require the same conditions (quiet conditions)!**

**If fluctuations in readings are larger than 15mmHg, and/or you hear irregular pulse tones, consult your doctor.**

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacture of your Microlife blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Esses (RWTüV-Essen).

**Never attempt to repair the instrument yourself!**

Any unauthorized opening of the instrument invalidates all warranty claims!

### Periodical recalibration

Sensitive measuring devices must from time to time be checked for accuracy. We therefore recommend a periodical inspection of the static pressure display every 2 years.

## 9. Warranty

---

Your Microlife blood pressure monitor is **guaranteed for five years** against manufacturers' defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

There are no user servicable parts inside. Batteries or damage from old batteries is not covered by the warranty.

**Please note: According to international standards, your monitor should be checked for accuracy every 2 years**

## 10. Reference to standards

---

### Device standard:

This device is manufactured to meet the European and United States standards for non-invasive blood pressure monitors:

EN1060-1 / 1995

EN1060-3 / 1997

EN1060-4 / 2004

AAMI/ANSI SP10

### Electromagnetic compatibility:

Device fulfills the stipulations of the European standard IEC60601-1-2

The stipulations of the EU-Guidelines 93/42/EWG for Medical Products Class IIa have been fulfilled.

## 11. Technical specifications:

---

Weight:	148g (with batteries + cuff)
Size:	85 x 77 x 75mm (including cuff)
Storage temperature:	-20 to +50°C (23°F - 122°F)
Humidity:	15 to 90% relative humidity maximum
Operation temperature:	10 to 40°C (50°F - 104°F)

Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	SYS/DIA: 30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Stores the last 30 measurements automatically for 2 Users
Measuring resolution:	1 mmHg
Accuracy:	Pressure within $\pm 3$ mmHg Pulse $\pm 5$ % of the reading
Power source:	2 dry cells (Batteries) size AAA 1.5V
Accessories:	storage case, PC software
Technical alterations reserved!	

## 12. How to contact us:

---

Distributed by:  
**Microlife USA, Inc.**  
 Dunedin, FL 34698

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