Blood Pressure Monitor

with Microlife Averaging Mode

Instruction Booklet for model #BP3AA1-2

Microlife Corporation is a proud national sponsor of the American Heart Association's High Blood Pressure Heart Profilers™

US Patent #'s 6,485,429 (Irregular Heartbeat Detector)
6,447,457 (Microlife Averaging Mode)

microlife®

Table of contents

1. Introduction

- 1.1 Features of your blood pressure monitor
- 1.2 Important information about self measurement

2. Important information on the subject of blood pressure and its measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?
- 2.4. Microlife Averaging Mode Technology

3. Components of your blood pressure monitor

4. Operation of your blood pressure monitor

- 4.1. Inserting the batteries
- 4.2. Using an AC power adapter
- 4.3. Cuff connection
- 4.4. Setting the time and date
- 4.5. Reading the set date

5. The measurement procedure

- 5.1. Before measuring
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measurement mode selection (single or average mode)
- 5.5. Setting memory for user 1 or user 2
- 5.6. Measuring procedure
- 5.7. Discontinuing a measurement
- 5.8. Memory – recall of the measurements
- 5.9. Clearing the memory

6. Error messages/Troubleshooting

7. Care and maintenance

8. Warranty

9. Certifications

10. Technical specifications

11. How to contact us

1. INTRODUCTION

1.1. Features of your Blood Pressure Monitor

Your Microlife blood pressure monitor is a fully automatic digital blood pressure measuring device for use on the upper arm. It enables very fast and reliable measurement of the systolic and diastolic blood pressure as well as the pulse by way of the oscillometric method. This device offers clinically proven accuracy and has been designed to be very user friendly.

Before using, please read this instruction manual carefully. Then store the manual in a safe place for later reference. Please contact your physician for further questions on the subject of blood pressure and its measurement.

ATTENTION!

1.2. Important information about self-measurement

Self-measurement means control, not diagnosis or treatment. Your values must always be discussed with your physician. You should never alter the dosages of any medication without direction from your physician.

The pulse display is not suitable for checking the frequency of heart pacemakers!!

In cases of irregular heartbeat, measurements made with this instrument should only be evaluated after consultation with your doctor.

Electromagnetic Interference: This device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens, etc.). These can lead to temporary impairment of the monitor's measuring accuracy.

2. IMPORTANT INFORMATION ON THE SUBJECT OF BLOOD PRESSURE AND ITS MEASUREMENT

2.1 How does high or low blood pressure arise?

The level of blood pressure is determined in the circulatory center of your brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels are altered. Blood vessel width is altered by fine muscles in the blood vessel walls. The level of arterial blood pressure changes periodically during heart activity. During the “blood ejection” (Systole), the value is highest (systolic blood pressure value). At the end of the heart’s “rest period” (Diastole) pressure is lowest (diastolic blood pressure value). Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2 Which values are normal?

Blood pressure is too high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, **while at rest**. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

If your systolic blood pressure values are between 140 mmHg and 160mmHg and/or the diastolic blood pressure values between 90 mmHg and 95mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician.**

Which values are normal?

The following standards for assessing high blood pressure (without regard to age) have been established by the National Institutes of Health JNC7.

Category	Systolic (mmHg)	Diastolic (mmHg)
Normal Self-Check	<120	<80
Pre-Hypertension	120 - 139	80 - 89
Stage 1 Hypertension	140 - 159	90 - 99
Stage 2 Hypertension	>160	>100

• Further Information

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension.” Consult your physician.
- Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.

2.3 What can be done if regular high or low values are obtained?

- a) Consult your physician.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constrictions caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your physician for information about the possible origins of your own increased blood pressure values.
- d) There are measures you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

1. Eating habits:

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

2. Previous illnesses

- Consistently follow all medical instructions for treating illness such as: diabetes (diabetes mellitus or sugar diabetes), fat metabolism disorder, and gout.

3. Habits

- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (coffee, tea and chocolate)

e) Your physical condition

- After a medical examination, and with your doctor's approval and direction; exercise.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your physician before beginning your exercise routine. You must receive advice regarding the type and extent of exercise that is appropriate for you.

2.4 Microlife Averaging Mode Technology


Microlife Averaging Mode (MAM) Technology is a new concept for optimum reliability in self-measurement of blood pressure. A single measurement is rarely your real result. An average of three readings is medically proven to provide the most accurate result. With the push of a single button, the Microlife Averaging Mode blood pressure monitor takes three consecutive measurements, and calculates the average.

3. COMPONENTS OF YOUR BLOOD PRESSURE MONITOR

The upper arm cuff included with your monitor is for arm circumferences of 9.5" - 13.25" (20-32cm).

Please Note: Do not force the cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. Larger cuffs are available.

AC Adapter and Cuff Connection
(AC Adapter sold separately)


4. OPERATION OF YOUR BLOOD PRESSURE MONITOR

4.1 Inserting the batteries

After you have unpacked your device, insert the batteries. The battery compartment is located on the back side of the device (see illustration).

- Remove cover from the bottom plate, as illustrated.
- Insert the batteries (4 x size AA). Always use AA long life batteries or alkaline 1.5v batteries. Do not use rechargeable batteries.
- If a battery warning appears in the display, the batteries are discharged and must be replaced. 

Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- If the blood pressure monitor is left unused for long periods, please remove the batteries from the device.


Functional check: Hold the ON/OFF button down to test all the display elements. When functioning correctly all segments must appear.

4.2. Using an AC adapter (optional accessory sold separately)

It is possible to operate this blood pressure instrument with an AC adapter (output 6 V DC / 600 mA with DIN plug).

- Insert the AC adapter plug into the socket at the left side of the instrument.
- Plug the AC adapter into a power socket. Test that power is available by pressing the ON/OFF button.

Note: No power is taken from the batteries while the AC adapter is connected to the instrument. If power is interrupted during the measurement, (e.g. by accidental removal of the AC adapter from the socket) the instrument must be reset by removing the plug from its socket and re-inserting the connectors. To purchase an AC adapter, please call Microlife 1-800-568-4147 or visit our website, www.microlifeusa.com.

4.3. Cuff connection

Insert the cuff tube into the opening provided on the left side of the instrument, as shown in the diagram.

4.4. Setting the time and date

This blood pressure monitor incorporates an integrated clock with date display. Time and date are recorded along with pressure and pulse. After new batteries have been inserted, the clock begins to run from the following setting: 2000-01-01 12:00 O'clock. You must then re-enter the date and current time. Please proceed as follows (Example: Entering 2002-06-15 Time 09:30 AM):

1. Press the User/Time button for at least 3 seconds. The display now indicates the set year, The last two characters are blinking.

5. Press the User/Time button again. The last two characters (day) are now blinking.

2. The correct year can be entered by pressing the MEMORY button. Once for each change. (Example: 1 x press)

6. The current day can now be entered by pressing the MEMORY button. (Example: 14 x presses advances the day from the 1st to 15th)

3. Press the User/Time button again. The display now switches to the current date, during which the first character (month) blinks. Note: Holding the button down speeds up the procedure.

7. Press the User/Time button again. The display now switches to the current time, during which the first character (Hour) blinks.


4. The current month can now be entered by pressing the MEMORY button. Example: pressing 5 x advances 5 months)

8. The corresponding hour can now be entered by pressing the MEMORY button. (Example: 9 x presses)

9. Press the User/Time button again. The last two characters (Minutes) now blink.



10. The minutes can now be entered by pressing the MEMORY button. (Example: 30 x presses advances the time from 0 to 30 minutes)

4.5. Reading the set date

After all settings have been made, click the User/Time button once. The date is briefly displayed and then the time. The input is now confirmed and the clock begins to run.

5. MEASUREMENT PROCEDURE

5.1. Before measuring:

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before the measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Compare measurements taken at the same time of day, since blood pressure changes during the course of the day.

5.2. Common sources of error:



Note: Comparable blood pressure measurements always require the same conditions! Conditions should always be quiet.

- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured! Each 15cm (6 inch) difference in height between your heart and the cuff results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please contact Microlife USA. Large cuff sizes are sold separately.
Note: Only use clinically approved Microlife Cuffs!
- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause or after your arm has been held up in order to allow the accumulated blood to flow away.

Please note, clinical studies have proven the Microlife Averaging Mode (MAM) system with a 15 second wait period between the second and third measurements. If you decide to retake your MAM measurement again, be sure to wait at least one minute before repeating a MAM measurement.

5.3. Fitting the cuff

- a) Pass the end of the cuff through the flat metal ring so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.)
- b) Place the cuff over the left upper arm so that the tube is closer to your lower arm. Align the cuff artery mark (arrow) over the arm artery than runs down the inner side of the arm.
- c) Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately $3/4$ " to 1" (2 to 3 cm) above the elbow and that the tube is closer to the inner side of the arm.
- d) Tighten the cuff by pulling the end and close the cuff by affixing the velcro. Make sure that the cuff Index end falls between the OK lines when the cuff is in position. If it doesn't, you are using the wrong cuff size.
- e) There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed.
- f) Lay your arm on a table (palm upward) so the cuff is at the same height as the heart. Make sure that the tube is not kinked.
- g) Remain seated quietly for five minutes before you begin the measurement.

Comment:

If it is not possible to fit the cuff to the left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

5.4. Measurement mode selection

Clinical studies demonstrate taking multiple blood pressure readings and calculating an “average” is more likely to determine your true blood pressure. Your Microlife Averaging Mode unit allows you to switch the unit to an Averaging Mode setting that automatically takes multiple readings!

Averaging Mode Slide Switch:

- If you would like to take an Averaging Mode measurement, move the slide switch to the right.
- Averaging Mode takes 3 measurements in succession and calculates the result. An “A” symbol in the display indicates that the unit is set to the Averaging Mode.
- There will be 15 seconds resting time in-between each measurement. The unit will count down from 15 seconds. A beep will sound at the last 5 seconds to remind you that the unit is about to take another measurement.

One is for a single measurement Three is Average Mode Selection


Single mode

- If you would like to take a single measurement, please move the slide switch to the left.

5.5. Setting memory for user 1 or user 2

Note: Before proceeding make sure the unit is turned on.

This advanced blood pressure monitor allows you to track blood pressure readings for 2 individuals independently.

- Before taking your measurement, make sure you set the unit for the intended user. (User 1, User 2)
- The unit is set to User 1. Click the User/Time button to change to User 2. The symbol indicates user 1 or user 2.
- Click the User/Time button again to change into User 1.

5.6. Measuring procedure

After the cuff has been appropriately positioned, the measurement can begin:






- Press the ON/OFF button. The display will show many icons followed by 3 beeps indicating "ready to measure."
- Press the START button. The pump begins to inflate the cuff. The rising pressure in the cuff is shown in the display of your monitor.
- After a suitable inflation pressure has been reached, the pump stops and the pressure gradually falls. The cuff pressures are displayed. If the inflation pressure is not sufficient, the monitor automatically re-inflates.
- When the instrument detects a pulse, the heart symbol in the display starts to flash and a beep is heard for every heartbeat.
- A long beep is sounded when the measurement is complete. The systolic, and diastolic blood pressure and pulse rate now appear in the display.
- The measurement reading remains on the display until you switch off the instrument. If no button is pressed for a period of 5 minutes the instrument switches itself off in order to preserve the batteries.
- When the unit is set to the Microlife Averaging Mode setting, 3 separate measurements will take place in succession after which your result is calculated. There is a 15 second resting time in-between each measurement. A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin.

If one of the measurements causes an error message, it will be repeated one more time. If any additional error occurs, the measurement will be discontinued and error code displayed.

Expanding the MAM measurement to see the three individual readings - This function allows you to view the three individual measurements used to calculate the MAM measurement. After taking a MAM measurement hold the memory button for 3 seconds until you hear a short beep. (Do not hold the button longer than 7 seconds or you will delete all the readings in the memory.) Let go of the button and watch the screen. It will automatically scroll through the three measurements used in the MAM measurement.


5.7. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g. the patient feels unwell), the ON/OFF button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

5.8. Memory – recall of the measurements (Unit should be turned on first)

This blood pressure monitor automatically stores each of the last 30 measurement values separately for 2 users. By pressing the MEMORY button, the last measurement as well as previous 29 measurements (MR29, MR28, ...MR1) can be displayed one after the other. Data obtained in Average Mode is indicated by an "A" symbol. If the memory is full, the next measurement will be added to memory and the oldest reading (MR 1) will be deleted.


5.9. Clearing the Memory – cancellation of all measurements

Attention! Before you delete all readings stored in the memory, make sure you won't need to refer to the readings at a later date. Keeping a written record is prudent and may provide additional information for your next physician visit. A log is included to record a history of your measurements.


In order to delete all stored readings, press the MEMORY button for at least 7 seconds. The display will show the symbol "CL" and 3 short beeps will sound to indicate deletion of stored readings. (Unit must be turned on.)

6. ERROR MESSAGES/TROUBLESHOOTING

If an error occurs during a measurement, a corresponding error code is displayed (Example: Error no. 1)

Error Number	Possible Cause(s)
Err 1	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. This situation can occur if the tube becomes unplugged after the systolic blood pressure has been measured. Check the tube connections and repeat the measurement. Further possible cause: No pulse could be detected
Err 2	Unnatural pressure impulses impair the result of the measurement. Cause: The arm was moved during the measurement. Remain still and repeat the measurement.
Err 3	Inflation of the cuff takes too long. The cuff is not correctly seated or the hose connection is not sealed. See section 5.3, fitting the cuff.
Err 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.

Other possible malfunctions troubleshooting

If problems occur when using the device, the following items should be checked.

Malfunction	Remedy
The display remains empty when the instrument is on. Batteries have been installed.	<ol style="list-style-type: none">1. Check the polarity of the batteries(+/-)2. If the display is erratic or unusual, remove the batteries and re-install new batteries.
The pressure does not rise although the pump is running.	<ol style="list-style-type: none">1. Check the connection of the cuff tube and connect properly if necessary.
The instrument frequently fails to measure the blood pressure values, or the values measured are too low or too high.	<ol style="list-style-type: none">1. Ensure the cuff fits correctly.2. Make sure the cuff is not too tight. Ensure tight clothing such as a rolled up shirt sleeve is not exerting pressure on the arm above the measuring point. Take off articles of clothing if necessary.3. Measure your blood pressure again in complete peace and quiet.
Every measurement produces a different value although the instrument functions normally and the values displayed are normal.	PLEASE NOTE: BLOOD PRESSURE FLUCTUATES CONTINUALLY SO SUCCESSIVE MEASUREMENTS WILL SHOW SOME VARIABILITY.
The blood pressure values measured differ from those measured by your physician.	<ol style="list-style-type: none">1. Record daily values, time and date, then consult your physician. PLEASE NOTE: INDIVIDUALS VISITING THEIR PHYSICIAN FREQUENTLY EXPERIENCE ANXIETY WHICH CAN RESULT IN A HIGHER READING THAN ONE OBTAINED AT HOME UNDER MORE RESTFUL CONDITIONS.

FURTHER INFORMATION

- The level of blood pressure is subject to fluctuations even with healthy people. Comparable measurements always require the same conditions (Quiet conditions)!
- For licensing, the device has been subjected to strict clinical tests and supervised by experienced Cardiac Specialists.
- This device is manufactured to meet European standards under the supervision of the Technical Monitoring Association (RWTüVEssen).
- If you have any questions regarding the use of this blood pressure monitor, please call toll free 1-800-568-4147 or visit our web site at www.microlifeusa.com. **Never attempt to repair the instrument yourself!**
- Any unauthorized opening of the instrument invalidates all warranties!

7. CARE AND MAINTENANCE, ACCURACY CHECK

a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.


b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of stress through twisting or buckling.

c) Handle the tube carefully. Avoid stress from pulling, pinching and sharp edges.

d) Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. Carefully remove spots on the cuff with a damp cloth and soap-suds. **The cuff should not be washed in a washing machine or in a dishwasher!**

e) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.

f) **Never open the device!** This can negatively impact the calibration (accuracy).

g) **Periodic accuracy check**

Sensitive measuring devices must be checked for accuracy from time to time. We recommend a periodical inspection of your unit by an authorized dealer every 2 years. Your factory authorized dealer would be pleased to provide more information regarding calibration. Call 1-800-568-4147 to find out about getting your blood pressure monitor calibrated.

8. WARRANTY

Your microlife blood pressure monitor is **guaranteed for 5 years** against manufacturers' defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by the warranty. The cuff is warranted for one year from date of purchase.

Please note: According to international standards, your monitor should be checked for accuracy every 2 years.

If fluctuations in readings are larger than 15 mmHg, and/or you hear irregular pulse tones, and/or the Irregular Heartbeat Detector icon appears frequently, consult your doctor.

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacture of your blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Essen (RWTÜV-Essen).

Never attempt to repair the instrument yourself!

Any unauthorized opening of the instrument invalidates all warranty claims!

9. CERTIFICATIONS

Device standard:

This device is manufactured to meet the European and United States standards for non-invasive blood pressure monitors:

EN1060-1

EN1060-3

EN1060-4

AAMI/ANSI SP10

Electromagnetic compatibility:

Device fulfills the stipulations of the International standard IEC60601-1-2

The B.H.S. (British Hypertension Society) has graded this product A/A or "Recommended for clinical and home use." This is the highest grading available for blood pressure monitors. Please see the B.H.S. web site at www.bhsoc.org.

*Product using the identical measurement algorithm was tested by unaffiliated researches using B.H.S. study protocol. Results on file and available by request.

11. TECHNICAL SPECIFICATIONS:

Weight:	650g (with batteries)
Size:	124 (W) x 205 (L) x 81 (H) mm
Storage temperature:	-20 to +50°C
Humidity:	15 to 90% relative humidity maximum
Operation temperature:	10 to 40°C
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	Oscillometric
Pressure sensor:	Capacitive
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Automatically stores the last 30 measurements separated for 2 users
Measuring resolution:	1 mmHg
Accuracy:	
Pressure	within ± 3 mmHg
Pulse	± 5 % of the reading
Power source:	DC 6V / 600mA
	a) 4 size AA batteries, 1.5V
	b) AC adapter / output : 6 V DC 600mA (Sold separately)
Special accessories:	Cuff type S102 (12" – 16") for arm circumference 31-41 cm (Sold separately) Cuff type S101 (9.5"-13.25") for arm circumference 22-32 cm
Reference to standards:	EN 1060-1 /-3, NIBP - requirements DIN 58130, NIBP - clinical investigations ANSI / AAMI SP10, NIBP - requirements

* Technical alterations reserved!

12. Contacts for Support

Distributed by:
Microlife USA, Inc.
424 Skinner Blvd. - Suite C
Dunedin, FL 34698

Toll Free: 1-800-568-4147

Email Questions to: custserv@microlifeusa.com
www.microlifeusa.com

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