Life Fitness™

BELIEVE IN YOU™

Premium Wrist Blood Pressure Monitor

Model # BP3AX1-4U
ID # 3AX1-4U-LF (344536)
Instruction Manual

PLEASE NOTE:
THIS MEDICAL INSTRUMENT MUST BE USED ACCORDING TO
INSTRUCTIONS TO ENSURE ACCURATE READINGS.

Questions? Preguntas?

Call Toll Free: 1-866-464-6184

Premium Wrist Blood Pressure Monitor

Instruction Manual

Model #BP3AX1-4U

ID#3AX1-4U-LF (344536)

Table of Contents

1. Introduction

- 1.1. Features of your blood pressure monitor
- 1.2. Important information about self-measurement

2. Important Information about Blood Pressure and its Measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. How can you treat high blood pressure?
- 2.4. Benefits of Life Fitness AveragingMode Technology™

3. The Components of Your Blood Pressure Monitor

4. Using Your Blood Pressure Monitor

- 4.1. Inserting the batteries
- 4.2. Initial setup
- 4.3. Setting alarm times
- 4.4. Sleep mode

5. Taking a Measurement

- 5.1. Before the measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measurement mode selection (LifeFitness Averaging Mode Technology™ feature)
- 5.5. Measurement procedure
- 5.6. Discontinuing a measurement

6. Memory Recall and User Setting

- 6.1. Viewing measurements saved in memory
- 6.2. How to use the individual memory feature
- 6.3. Clearing the memory

7. Software Functions

- 7.1. Installation of software and downloading readings

8. Adjusting the Display Contrast

9. Error Messages / Malfunctions

10. Care and Maintenance

11. Warranty

12. What U.S. and International Standards is this device manufactured to meet?

13. Technical Specifications

14. How to Contact Us

1. Introduction

1.1. Features of your Premium Blood Pressure Monitor, ID # 3AX1-4U-LF

The blood pressure monitor is a fully automatic, digital blood pressure measuring device for use on the wrist. This enables very fast and reliable measurement of systolic and diastolic blood pressures as well as the pulse through the oscillometric measuring method. The device offers excellent accuracy proven in clinical studies and has been designed for ease of use.

Before using, please read this instruction manual carefully and then keep it in a safe place. For further questions on the subject of blood pressure and its measurement, please contact your doctor.

Caution!

1.2. Important information about self-measurement

- Self-measurement means knowledge, not diagnosis or treatment. Unusual values must always be discussed with your doctor. **You should never alter the dosages of any drugs prescribed by your doctor.**
- The pulse reading is not suitable for checking the frequency of heart pacemakers!
- Clinical trials have shown that wrist measurements can reach similar accuracy to upper arm measurements. However, some individuals may get a different result between both methods. Therefore, we generally recommend to compare the readings of this device to data obtained by your doctor.
- Make sure you keep the device at heart level during measurement to avoid a user accuracy error.

Electromagnetic interference

The device contains sensitive electronic components. Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave cookers). These can lead to temporary impairment of the measuring accuracy.

2. Important information about blood pressure and its measurement

2.1 How does high or low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of your brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is affected by fine muscles in the blood vessel walls. Your level of arterial blood pressure changes periodically during heart activity. During the "blood ejection" (Systole), the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value). Your blood pressure values should lie within certain normal ranges in order to prevent particular diseases.

2.2 Which values are normal?

Blood pressure is dangerously high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, while at rest. In this case, please consult your physician immediately. Long-term values at this level endanger your health through continual damage to the blood vessels in you body. If your systolic blood pressure values lie between 140 mmHg and 160mmHg and/or your diastolic blood pressure values lie between 90 mmHg and 95mmHg consult your physician. Regular self-checks will be necessary. If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician. Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately. If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician.**

The following standards for assessing high blood pressure (without regard to age) have been established by the U.S. National Institutes of Health JNC7, 2003.

Category	Systolic (mmHg)	Diastolic (mmHg)
Normal Self-Check	<120	<80
Pre-Hypertension	120 - 139	80 - 89
Stage 1 Hypertension	140 - 159	90 - 99
Stage 2 Hypertension	>160	>100

Further information

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called "labile hypertension." Consult your physician.
- Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.

2.3 What can be done if regular high or low values are obtained?

- a) Consult your physician.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body may be endangered due to narrowing caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between common primary (essential) hypertension, and secondary hypertension. Secondary hypertension can be ascribed to specific organ malfunctions. Please consult your physician for information about the possible origins of your own high blood pressure.
- d) There are measures you can take to reduce and even prevent high blood pressure. These measures should include permanent lifestyle changes.

1) Eating habits:

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of salt.
- Avoid fatty and fried foods.

2) Previous illnesses

- Consistently follow all medical instructions for treating illness such as: diabetes (diabetes mellitus or sugar diabetes), fat metabolism disorder, and gout.

3) Habits

- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea and chocolate).

4) Your physical condition


- After a medical examination, and with your doctor's approval and direction, exercise.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your physician before beginning your exercise routine. You must receive advice regarding the type and extent of exercise that is appropriate for you.

2.4 Life Fitness AveragingMode Technology™

Life Fitness AveragingMode Technology™ is a new concept for optimum reliability in measurement of blood pressure. A single measurement is rarely your real result. An average of three readings is medically demonstrated to provide the most accurate result. With the push of a single button, the Life Fitness AveragingMode Technology™ blood pressure monitor takes three consecutive measurements and calculates the average.

3. Components of your blood pressure monitor

The illustration below shows your blood pressure monitor's key features:


This wrist cuff fits wrist circumference of 5.3" to 7.67" (13.5cm to 19.5 cm).

4. Using the blood pressure monitor

4.1. Inserting the batteries

After you have unpacked your device, insert the batteries.

1. Remove the cover as illustrated (A).
2. Lay the clear plastic strip (B) inside the battery compartment, and insert 2 x size AAA 1.5V batteries (C), observing the indicated polarity according to the markings on the cuff (D).

Caution!

- After the low battery warning appears, the device is disabled until the batteries have been replaced.
- If the battery warning appears in the display, the batteries are discharged and must be replaced by new ones.
- Use AAA Long-Life or Alkaline 1.5V Batteries. The use of rechargeable batteries is not recommended.
- If the blood pressure monitor is left unused for long periods, remove the batteries from the device.

4.2. Initial Setup

When the batteries are installed, the display shows the language selection screen.

Press **▲** or **▼** to select the display language. Press **OK (M)** to confirm the selection. The default language is English.

The device prompts you to set the time and date. The year blinks highlighting the first field to change.

Press **▲** or **▼** to change the year. Press **OK (M)** to confirm the date. The month, day, hour, and minute are all consecutively highlighted. Press **▲** or **▼** to change these fields and press **OK (M)** to confirm. You are then prompted to set the alarms.

4.3. Setting Alarm Times

After you have set the time, you can set one or two alarms to alert you to take your blood pressure measurements. An alarm icon appears in the display. Press **▲** or **▼** to set the hour and minute for Alarm 1. Press **OK (M)** to confirm each setting.

When the alarm icon is highlighted, press **▲** or **▼** to enable  or disable  Alarm 1. Press **OK** to confirm. You are then prompted to set Alarm 2. Follow the instructions for Alarm 1 to set this alarm.

After setting the alarms, press **OK** to exit the Time and Date Setup. An alarm icon  appears in the display to indicate that the alarm is activated.

Attention!


The time display is set on a 24-hour clock.

When the final field (DISPLAY OFF) is highlighted, press ▲ or ▼ to select whether you want the time and date display to remain on when the device is switched off. Press **OK** to confirm the settings and exit the time setup screen.

4.4. Sleep Mode

The device is designed to stay in sleep mode as batteries are loaded. The illustration shows the display elements, when the BPM is in sleep mode, and the DISPLAY OFF feature is disabled (see Section 4.2).

5. Taking a measurement

5.1. Before the measurement

- Avoid eating, smoking as well as all forms of exertion directly before the measurement. All these factors influence the measurement result. Try to find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement on the left wrist.
- If you compare readings note the time of day as blood pressure fluctuates throughout the day.

5.2. Common sources of error:

Note: Comparable blood-pressure measurements always require the same conditions! Quiet conditions are always required for your optimal measurement.

- All efforts by the patient to support the arm can increase the blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion

for support if necessary.

- If the wrist artery lies considerably lower or higher than the heart, an erroneously higher or lower blood pressure will be measured! (Each 15cm or 5.9" difference in height results in a measurement error of 10mmHg!)
- A loose cuff causes false measurement values.
- With repeated measurements, blood can accumulate in the respective arm, which can lead to false results. Correctly executed blood pressure measurements should therefore first be repeated after a 1 minute pause or after the arm has been held in order to allow the accumulated blood to flow away. (Does not apply to a Life Fitness AveragingMode™ measurement.)

5.3. Fitting the cuff

- a) Remove all objects and jewelry (e.g. wristwatch) from the wrist in use. Put the cuff over the wrist with the monitor face on the inside of the wrist.
- b) The distance between the cuff and the hand should be approx. 1/2" (10 mm).
- c) Secure the cuff with the Velcro™ fastener, so that it lies comfortably and not too tight. No space should remain between the cuff and the wrist.
- d) Lay your arm on a table, with your palm upwards. Support your arm with a cushion, so that the cuff rests at about the same height as your heart. Do not hold or support the cuff with your free hand. Remain still for 5 minutes before beginning a measurement.


5.4. Measurement Mode Selection

Choose which measurement mode you want to use before taking a blood pressure reading. Press **OK (M)** to enter the menus, press **▲** or **▼** to select Measuring mode, and press **OK (M)** to enter the submenu Select , «1» or «3/MAM» for Life Fitness AveragingMode Technology™ with the scroll buttons and confirm your selection by **«OK»**.

Life Fitness AveragingMode Technology™

In average mode, the unit takes 3 continuous measurements, and afterwards calculates the result based on an average using Life Fitness AveragingMode Technology™ (see Section 2.4). There is a 15-second rest time between measurements; the unit beeps during the last 5 seconds to remind you of the upcoming measurement.

Single Mode

In single mode «1», only one measurement is taken.


5.5. Measurement Procedure

After the cuff and arm have been appropriately positioned and the measurement mode is set, the measurement can be taken as follows. Press the ON/OFF/START button. The Life Fitness Welcome screen appears, and is followed by a screen which shows you how to properly hold your arms. If the display off function is enabled, the Sleep Mode screen appears for two seconds before measurements are taken.

The pump begins to inflate the cuff and the increasing cuff pressure is displayed.


Hold left wrist at heart level and remain still during measurement.


After the suitable inflation pressure has been reached, the pump stops, the measurement signals are displayed and the heart icon flashes.

This curve shows each heart beat activity. Movement artifacts (excessive movement) would become visible by large deviations within the curve.


A long beep sounds when the measurement has been completed. The systolic and diastolic blood pressure and pulse rate appear in the display. The measurement reading remains on the display for approximately 30 seconds. If no button is pressed during the 30-second period the instrument switches itself off in order to preserve the batteries.

When the unit is set to MAM, it means the Life Fitness AveragingMode Technology™ has been engaged and three separate measurements take place as explained under section 5.4.

Countdown timer

Measurement number

5.6. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason, for example, if the patient feels ill, press the ON/OFF/START button. The device immediately lowers the cuff pressure and enters into sleep mode.

6. Memory handling and user setting

6.1 Viewing Measurements Saved in Memory

This blood pressure monitor automatically stores each of the last 30 measurement values. To view the stored measurements for a user, press ▲ or ▼ in sleep mode. The last measurement taken for the user is displayed.

The user, measurement sequence number, date, and time are displayed on the top.

Data obtained in Life Fitness AveragingMode Technology™ are indicated by a MAM symbol in the upper right corner.

- Press ▲ or ▼ to scroll through the list items.
- The first item in the list is the average of all stored data.
- Press **OK** and the display shows table of all data.
- The first entry is the average. Data obtained in Life Fitness AveragingMode Technology™ is indicated by an «A» symbol at the very right.
- Press ▲ or ▼ to scroll through the list.


Life Fitness
AveragingMode
Technology™

Measurement
number

	9	MAM
3/06		15:48
sys	124	mmHg
dia	73	mmhg
pul	77	/min.

	3/06	15:48	
No.	sys	dia	pul
avg.	149	110	85
1	145	109	83
2	119	71	74 A
3	124	73	79
4	130	80	77 A

Press **OK** and the display shows graph of all data. The lines indicate the Systolic (top end) and Diastolic (bottom end) values. Press **▲** or **▼** to move the cursor between measurements. The currently selected measurement is indicated by a blinking line. Press **ON/OFF/START** to exit this mode.

6.2. Setting the User

This monitor allows you to track readings for two individuals independently. Scroll to User and press **OK** to enter the submenu. Press **▲** or **▼** to make your selection. Select

 to set the device to not save readings.

Select  to select user 1. Any readings made when user 1 is selected will be stored only in memory 1.


Select  when user 2 is selected will be stored only in memory 2.

Note: There is a User 0, which does not store memories. If your unit is on the User 0 setting, no readings will be stored in memory.


After selecting user 1 or user 2 you are prompted to enter the user name. You can use up to six characters.

To keep the default name «User_1», press **ON/OFF/START**. To enter a different name, follow these instructions:

Press **▲** or **▼** to scroll through the alphanumeric characters A ~ Z and 0 ~ 9. There is also one underscore and one space

*A friend can measure their blood pressure without saving it in memory by setting the user to the  symbol (User 0).

character.

Press **OK** (M) to accept the character and move to the next space.

After you have finished entering the name, press **OK** to exit the user setup menu. The user name is displayed on the screen.

6.3. Clearing the Memory



This item enables you to clear readings for user 1 or user 2. You can erase single readings, or erase all readings at once. Select user 1 or user 2 as explained under section 6.2. Scroll to Clear Memory and press **OK**. The user you selected is shown in the upper left corner.

Clear Single Readings

Select «SINGLE» by the selection buttons and press **OK** to erase single readings. Press ▲ or ▼ to select a reading from the list. An arrow ► highlights the selected reading. Press **OK** to erase the selected reading.


Clear All Readings

Select «All» by the selection buttons and press **OK** to erase all readings of the selected user. You are prompted to confirm the action. Press **OK** to erase all readings of the selected user. Press the ON/OFF/START button to cancel deletion of all readings.

CLEAR MEMORY
3/06 15:48

No.	SYS	DIA	PUL
1	124	84	76
2	139	95	74
3	139	98	80
4	146	106	72

7. Software Functions

(For more detailed information, see instructions included on the CD)

Measurements stored in your wrist unit can be downloaded to your PC using the enclosed software. Your PC will allow a capacity of monitoring 80 patients, each with 1000 data (note: overuse will lower system efficiency). The memory data can be transferred to the PC by connecting the monitor via the USB cable. **Important: your wrist blood pressure monitor must be set to PC LINK mode whenever the unit is in contact with your computer.** When the software prompts you to make the connection to your blood pressure monitor, your monitor must be in the PC LINK mode. Push the M Button and choose PC LINK then push the M Button again to select it. You must set your wrist blood pressure monitor to PC LINK mode to install the software and every time you want to download your measurements to your PC. *When the unit is in PC LINK mode you will see a series of curved lines that pulsate from a stationary dot. Note: Software does not work with MAC computers!*

System Requirements for Blood Pressure Analyzer Software:

- Windows 98 SE, 2000, XP
- CD-Rom Drive
- Minimum 256MB RAM
- 500MB Available Hard Disk Space
- USB Port version 1.0 or higher
- Microlife Blood Pressure Monitor with USB Port

7.1. Installation and Data Transmission

- Insert CD into CD ROM drive of your PC. The installation will start automatically (If not, click on SETUP.EXE). The software will prompt you to connect your unit to the computer via the USB cable. The wrist unit must be in PC LINK mode.
- To transmit data you must have stored data in the memory first. To download memories from User 1, you must set the blood pressure monitor to User 1 prior to downloading. To download memories from User 2, you must set the monitor to User 2. There is a User 0, which does not store memories. If your unit is on the User 0 setting, no readings will be downloaded. Remember to make sure your wrist unit is in PC LINK mode prior to connecting to your computer.
- You cannot take a blood pressure measurement while the unit is in PC LINK mode.

During the PC connection, the device is completely controlled by the computer. Please refer to the 'help' file for software instructions.


Note: For problems with installation, go to Error Messages/Troubleshooting in

Section 9.

8. Adjusting the Display Contrast

This item enables you to adjust the contrast. Scroll to Contrast and press **OK** to enter the submenu. Press  to increase the display contrast. Press  to decrease the display contrast. After you have adjusted the contrast, press **OK** to confirm the setting.

Note: Do not set the contrast to the lowest setting. It will appear blank and remain that way next time you turn it on, rendering it impossible to see. Make only minor contrast adjustments.

9. Error messages / malfunctions

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed (Example: Error No. 2).

Error No.	Possible cause(s)
ERR 1	No pulse was detected.
ERR 2	Unnatural pressure impulses influence the measurement result. Reason: Your arm was moved during the measurement (Referred to as an "artifact").
ERR 3	Inflation of the cuff takes too long. The cuff is not correctly seated.
ERR 5	The measured readings indicated an unacceptable difference between systolic and diastolic pressures. Take another reading following directions carefully. Contact your doctor if you continue to get unusual readings.
ERR 6	Single data differ too much during Life Fitness AveragingMode™ (MAM) even after 3 cycles. No average result can be displayed.

Other possible malfunctions and their elimination

If problems occur when using the device, the following points should be checked and if necessary, the corresponding measures are to be taken:

Malfunction	Remedy
The display remains blank when the instrument is switched on although the batteries are in place.	<ol style="list-style-type: none">1. Check batteries for correct polarity and if necessary insert correctly.2. If the display is unusual, re-insert batteries or exchange them.3. Be sure contrast is not turned too far down.
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none">1. Check the positioning of the cuff.2. Measure the blood pressure again in peace and quiet. See section 5, Taking a Measurement
Every measurement produces a different value although the instrument functions normally and the values displayed are normal	<ol style="list-style-type: none">1. Please read the following information and the points listed under «Common sources of error». Repeat the measurement. <p>Please note: Blood pressure fluctuates continually so successive measurements will show some variability.</p>
Blood pressure measured differs from those values measured by the doctor.	<ol style="list-style-type: none">1. Record the daily development of the values and consult your doctor. <p>Please note: Individuals visiting their doctor frequently experience anxiety which can result in a higher reading at the doctor's office than at home under resting conditions.</p>
Trouble downloading measurements to my PC	Uninstall the software using uninstall feature from the cd. Take the cd out. Restart PC. Reinstall software. Turn on your wrist monitor and put into PC Link mode. Make sure you have measurements stored in User 1 or User 2 then set on appropriate user (Not User 0). Attach wrist monitor to computer via supplied USB cable and click OK on computer screen.

10. Care and maintenance

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle the cuff carefully and avoid all types of stress through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gasoline, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. The cuff must not be washed!
- d) Do not drop the instrument or treat it roughly in any way. Avoid strong vibrations.
- e) Never open the device! Otherwise the manufacturer's warranty becomes invalid!

Further information

Blood pressure is subject to fluctuations even in healthy people.

Comparable measurements always require the same conditions (quiet conditions)! If fluctuations in readings are larger than 15mmHg consult your doctor.

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacturing of your Life Fitness blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Esses (RWTüV-Essen).

Never attempt to repair the instrument yourself!

Any unauthorized opening of the instrument invalidates all warranty claims.

Periodical accuracy verification

Sensitive measuring devices must from time to time be checked for accuracy.

We therefore recommend a periodical inspection of the static pressure display every 2 years.

11. Warranty

Your Life Fitness blood pressure monitor is **guaranteed for a lifetime** against manufacturer's defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

There are no user servicable parts inside. Batteries or damage from old batteries is not covered by the warranty.

Please note: According to international standards, your monitor should be checked for accuracy calibration every 2 years

12. Reference to standards

Device standard:

Device corresponds to the requirements of the European standard for non-invasive blood-pressure monitor
EN1060-1 / 1995
EN1060-3 / 1997
EN1060-4 / 2004

AAMI/ANSI SP10

Electromagn. compatibility:

Device fulfills the stipulations of the European standard IEC60601-1-2

The stipulations of the EU-Guidelines 93/42/EWG for Medical Products Class IIa have been fulfilled.

12. Technical specifications

Weight:	140g (with batteries)
Size:	67 x 70 x 80 mm
Storage temperature:	-20 to +50°C (23°F - 122°F)
Humidity:	15 to 90% relative humidity maximum
Operation temperature:	10 to 40°C (50°F - 104°F)
Display:	dot-matrix LCD (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	solid-state
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 per minute
Measuring resolution:	1 mmHg
Accuracy:	Pressure within ± 3 mmHg Pulse ± 5 % of the reading
Memory:	Storing the last 30 measurements automatically for each of two users
Power source:	2 dry cells (Batteries) UM-4, size AAA 1.5V
Accessories:	Storage case

Technical and color alterations may differ.

13. How to Contact Us

Advanced Healthcare Distributors, LLC
Concord, NH 03301

Toll Free: 1-866-464-6184

Email Questions to: LifeFitnessSupport@microlifeusa.com

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