Premium Portable Blood Pressure Monitor

with Irregular Heartbeat Detection,
and Microlife Averaging Mode

Instruction Booklet for model# BP3MQ1

US Patent #'s 6,485,429 (Irregular Heartbeat Detector)
6,447,457 (Microlife Averaging Mode)

American Heart
Association®

Learn and Live...

Microlife Corporation is a proud national
sponsor of the American Heart Association's
High Blood Pressure Heart Profilers™

microlife®

Automatic Premium Portable Blood Pressure Monitor - Instruction Manual

Table of contents

1. Introduction

- 1.1. Features of your blood pressure monitor
- 1.2. Important information about self measurement

2. Important information on blood pressure and its measurement

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done if regular high/low values are obtained?

3. Components of your blood pressure monitor

4. Using your monitor for the first time

- 4.1. Inserting the batteries
- 4.2. Using an AC Adapter
- 4.3. Setting the date and time
- 4.4. Cuff tube connection
- 4.5. Select the user
- 4.6. Setting the medication reminder
- 4.7. Microlife Averaging Mode

5. Measurement Procedure

- 5.1. Before measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Irregular Heartbeat Detector
- 5.6. Memory - Displaying the last 99 measurements
- 5.7. Discontinuing a measurement
- 5.8. Battery Charge Indicator

6. Error Messages/Troubleshooting

7. Care and maintenance

- 7.1. Accuracy Test

8. Warranty

9. Certifications

10. Technical specifications

11. How to contact us

1. Introduction

1.1. Features of your Microlife Premium Portable Automatic Blood Pressure Monitor

Your Microlife blood pressure monitor with Microlife Averaging Mode and Irregular Heartbeat Detection is a fully automatic digital blood pressure measuring device for use by adults on the upper arm at home or in your doctor's/nurse's office. It enables very fast and reliable measurement of systolic and diastolic blood pressure as well as pulse by way of the oscillometric method. This device offers clinically proven accuracy and has been designed to be portable and very user friendly.

Before using, please read this instruction manual carefully and then keep it in a safe place. Please contact your doctor for further questions on the subject of blood pressure and its measurement.

Attention!

1.2. Important information about self-measurement

- **Self-measurement means control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.**
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In cases of irregular heartbeat (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with your doctor.

Electromagnetic interference:

The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy.

2. Important information on blood pressure and its measurement

2.1. How does high/low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is affected by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the "blood ejection" (Systole) the value is highest (systolic blood pressure value). At the end of the heart's "rest period" (Diastole) pressure is lowest (diastolic blood pressure value). Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

2.2. Which values are normal?

Blood pressure is very high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, **while at rest**. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in you body.

If your systolic blood pressure values are between 140 mmHg and 159mmHg and/or the diastolic blood pressure values between 90 mmHg and 99mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician.**

Which values are normal?

The following standards for assessing high blood pressure (without regard to age) have been established by the National Institutes of Health JNCV7, 2003.

Category	Systolic (mmHg)	Diastolic (mmHg)
Normal	<120	<80
Pre-Hypertension	120-139	80-89
Hypertension		
Stage 1 Hypertension	140-159	90-99
Stage 2 Hypertension	>160	>100

Further information

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension.” Consult your doctor.
- **Correctly measured diastolic blood pressure values above 120mmHg require immediate medical treatment.**

2.3. What can be done if regular high or low values are obtained?

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

A) Eating habits

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

B) Previous illnesses

- Consistently follow all medical instructions for treating illness such as:
 - Diabetes (Diabetes mellitus or sugar diabetes)
 - Fat metabolism disorder
 - Gout

C) Habits


- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea, chocolate).

D) Your physical condition

- **After a medical examination, and with your doctor's approval and direction, exercise.**
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. You must receive advice regarding the type and extent of exercise that is appropriate for you.

3. Components of your portable blood pressure monitor

a) Measuring unit


*AC Adapter sold separately. Call 1-800-568-4147 to purchase AC Adapter Model #HMA-AC2

b) Upper arm cuffs:

Type S101 for arm circumference 24–33.7 cm or 9.5" - 13.25"

Type S102 for arm circumference 30.5 - 40.64 cm 12" - 16"

Your monitor comes with a medium and a large sized cuff.

Please Note:

Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call 1-800-568-4147 for further information. You may prefer to use one of our wrist type monitors.

4. Using your monitor for the first time


4.1 Inserting the batteries

After you have unpacked your device, insert the batteries. The battery compartment is located on the back of your device. (See illustration)

- Remove cover as illustrated.
- Insert the 4 AA 1.5V batteries following the indicated polarity.
- If the battery warning appears on the display, the batteries are discharged and must be replaced.

Attention!

- After the low battery indication warning appears, the device will not function until the batteries have been replaced. Please use "AA" Long Life or Alkaline 1.5V batteries. Rechargeable batteries are not recommended.
- If the blood pressure monitor is left unused for long periods of time, remove the batteries from the device.

4.2 Using an AC Adapter

(Special accessory sold separately. This unit requires a special AC Adapter. Call 1-800-568-4147 to purchase AC Adapter Model #HMA-AC2)

It is possible to operate this blood pressure monitor with an AC adapter. Once you receive your AC Adapter in the mail, make certain that you use an adapter that fulfills the legal requirements and electronic requirements in the U.S. (UL standard).

- a) Push the plug into the back side of the instrument.
- b) Push the AC adapter into a 110V power socket (US or Canada).

Note:

- No power is taken from the batteries while the AC adapter is connected to the device.
- If the power is interrupted during the measurement (e.g. by removal of the adapter from the wall socket), the instrument must be reset by removing the plug from the device.
- Please call 1-800-568-4147 if you have questions relating to the AC Adapter (Model # HMA-AC2).

4.3. Setting the date and time

With the unit OFF, press the Time/Date button on the side of the monitor and hold for three seconds until you see the  icon start flashing in the upper left corner of the monitor screen. The Memory button becomes your scroll button and the Time/Date button moves you to the next time, date or user.

When the User icon is flashing, you can toggle between User 1 and User 2 by pressing the M (Memory button). If you press the M (Memory button) one time and the User changes to User 2, you can press the Time/Date button to set your unit on User 2 or you can press the M (Memory button) again to go back to User 1. Either way you must press the Time/Date button on the side of the monitor to confirm the setting and move on. The unit will automatically move on to the Time/Date settings. If you don't want to change the Time and Date then just keep pressing the Time/Date button six times until nothing is flashing on the monitor screen.

If you want to change the Date and Time, then press the Time/Date button once the correct User is flashing on the screen so you can move on to set the correct year.

The unit will automatically move to the Year setting after the User has been set. To change the year, press the M (Memory button) to scroll forward. **Please note:** You cannot

Time/Date button
on left side of
monitor.

scroll back so scroll slowly so you don't miss the correct year. If you do pass the correct year the unit will scroll all the way to 2050, then start again at the year 2000. Once you select the correct year, press the Time/Date button to set it and move on to set the correct month. Then the Month will begin flashing. Again, use the M (Memory button) to select the correct month and press the Time/Date button to set it and move on to the day setting.

After you set the month, the Day will begin flashing. Again, use the M (Memory button) to select the correct day and press the Time/Date button to set it and move on to the hour setting.

After you set the day, the Time hour will begin flashing. Again, use the M (Memory button) to select the correct hour according to AM or PM and press the Time/Date button to set it and move on to the minutes.

The last setting is the Time Minutes. The minutes will be flashing. Again, use the M (Memory button) to select the correct minutes and press the Time/Date button to set it and complete your setting.

Please NOTE: You can cancel out of setting the User/Time/Date by pressing the Start/Stop button. If no operation occurs during the User/Time/Date setting process, the unit will automatically cancel out of it. To re-initiate the setting process, press and hold the Time/Date button for three seconds until the User icon begins flashing. Start the procedure all over again or press the Time/Date button again and again until you reach the desired location to reset.

4.4. Cuff tube connection

Insert the cuff tube into the opening on the left side of the monitor indicated by a drawing of the cuff (see drawing at right).

4.5. Select the User

This blood pressure monitor is designed to store 99 measurements for each of two users. Before taking a measurement, be certain that the correct user has been selected.

- a) With the unit off, press and hold the Time/Date button for 3 seconds until the user icon in the upper left corner of the LCD screen flashes.
- b) Press the Memory button to toggle between users.
- c) Press the START/STOP button to make your selection.

4.6. Setting the medication reminder


This instrument allows you to set two alarm times at which an alarm signal will then be triggered. This can be a useful aid, for instance as a reminder to take medication or to remind you to take your blood pressure at the same time each day.

1. To set an alarm time, press the Time/Date button on the side of the unit (the instrument must have been switched off beforehand) and immediately afterwards the "M" (Memory button) and hold both down until the bell symbol appears in the bottom left of the display. Then release both buttons. The flashing "1" in the display indicates that the first alarm time can now be set.
2. Press the Time/Date button on the side of the unit to set the hours-the hours display flashes and pressing the "M" (Memory button) allows you to set the alarm hour. To confirm, press the Time/Date button.
3. The minute display will now flash. The minutes can be set using the "M" (Memory button). To confirm, press the time button again.
4. The bell symbol will now flash. Use the "M" (Memory button) to select when the alarm time is to be active (Bell icon) or inactive (when inactive a Bell icon with a cross over it appears). To confirm, press the Time/Date button on the side of the unit.

To set a second alarm (Medication or Blood Pressure Measurement Reminder), proceed as above but if the "1" flashes, press the "M" button to select "2" and confirm with the Time/Date button on the side of the unit.

- An active alarm time is indicated by the bell symbol in the display.
- The alarm sounds at the selected time each day.
- To switch-off the alarm when it is sounding, press the Start/Stop button.
- To permanently switch-off the alarm, proceed as above (steps 1-4) and select the crossed-out bell symbol. This will then disappear from the display.
- The alarm times must be re-entered each time the batteries are replaced.


4.7. Microlife Averaging Mode

In Microlife Averaging Mode, three measurements are automatically taken in succession and the result then automatically analyzed and displayed. Because your blood pressure constantly fluctuates, a result determined in this way is more reliable than one produced by a single measurement.

The switch on the side of the blood pressure monitor (right above the Time/Date button but below the AC Adapter connection port [AC Adapter sold separately]) allows you to select whether you want to take a single measurement or three consecutive measurements in Microlife Averaging Mode. If you slide the switch toward the 3, the unit is in averaging mode. The resulting measurement on the display will indicate that the measurement is from Microlife Averaging Mode because an icon  will appear when a measurement was taken in this mode.

After pressing the START/STOP button the selected Microlife Averaging Mode appears in the display as the  symbol.

Single Mode: If you want a single measurement then slide the switch toward the 1 on the side of the unit. The single measurement will take just a single reading.

*AC Adapter sold separately. Call 1-800-568-4147 or visit www.microlifeusa.com to purchase the special AC Adapter (Model #HMA-AC2) for this unit.

5. Measurement Procedure

Please note: You should always be seated before and during measurement.

5.1. Before measurement:

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always compare measurements taken at the same time of day, since blood pressure changes during the course of the day.

5.2. Common sources of error:

Note: Comparable blood pressure measurements always require the same conditions! Conditions should always be quiet.



- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured! Each 15cm (6 inch) difference in height between your heart and the cuff results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please call 1-800-568-4147.

Note: Only use clinically approved cuffs!

- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause or after your arm has been held up in order to allow the accumulated blood to flow away.

5.3. Fitting the cuff

- Pass the end of the cuff through the flat metal ring so that a loop is formed. The Velcro™ closer must be facing outward. (Ignore this step if the cuff has already been prepared.)
- Place the cuff over the left upper arm so that the tube is closer to your lower arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 3/4" to 1" (2 to 3cm) above the elbow and that the tube is closer to the inner side of the arm.


Important! The **small white arrow (Artery Mark)** on the cuff must lie exactly over the artery which runs down the inner side of the arm.

- Tighten the cuff by pulling the end and close the cuff by affixing the Velcro™.
- There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed. Cuffs that don't fit properly result in false measurement values. Measure your arm circumference if you are not sure of proper fit. Larger cuffs are available (Section 3B).
- Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.

ARTERY MARK

- Remain seated quietly for at least two minutes before you begin the measurement.

Comment:

If it is not possible to fit the cuff to your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

Comparable blood pressure measurements always require the same conditions (Relax for several minutes before a reading).

5.4. Measuring procedure

After the cuff has been appropriately positioned on the arm and connected to the blood pressure monitor the measurement can begin:

- Press the START/STOP button. The pump begins to inflate the cuff. In the display, the increasing cuff pressure is continually displayed.
- After automatically reaching an individual pressure, the pump stops and the pressure slowly falls. The cuff pressure is displayed during the measurement.
- When the device has detected your pulse, the heart symbol in the display begins to blink.
- When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values, as well as the pulse are now displayed.
- The appearance of this symbol  signifies that an irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements. NOTE: We recommend contacting your physician if you see this indicator frequently.
- The measurement results are displayed until you switch the device off. If no button is pressed for 1 minute, the device switches off automatically.
- You can take a measurement while holding this monitor in your hand if you do not want or have a hard surface to place the unit on. However, you should always hold the unit in the opposite arm of the cuff. If the cuff is on the left arm, hold the blood pressure monitor in your right hand.

START/STOP
Button

Pumping
Pressure


Measuring

Systolic

Diastolic

Pulse

Measurement
complete

When the unit is set to the Averaging Mode setting, 3 separate measurements will take place in succession after which your result is calculated. There is a 15 second resting time in-between each measurement. A count down indicates the remaining time and a beep will sound 5 seconds before the 2nd and 3rd readings will begin.


If one of the measurements causes an error message, it will be repeated one more time. If any additional error occurs, the measurements will be discontinued and an error code displayed.

Expanding the averaged measurement to see the three individual readings:

This function allows you to view the three individual measurements used to calculate the Microlife Averaging Mode (MAM) measurement. After taking a MAM measurement hold the memory button for 3 seconds until you hear a short beep. (Do not hold the button longer than 7 seconds or you will delete all the readings in the memory.) Let go of the button and watch the screen. It will automatically scroll through the three measurements used in the measurement.

5.5. Irregular Heartbeat Detector

This symbol  indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal basal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times a week with measurements taken daily), we advise you to tell your doctor. Please show your doctor the following explanation:

Information for the doctor on frequent appearance of the Irregular Heartbeat Symbol

This instrument is an oscillometric blood pressure monitor device that also analyzes pulse frequency during measurement. The instrument is clinically tested.

If pulse irregularities occur during measurement, the irregular heartbeat symbol is displayed after the measurement. If the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

5.6. Memory – displaying the last 99 measurements

At the end of a measurement, this monitor automatically stores each result with date and time. This unit stores 99 memories for each of 2 users.

Viewing the stored values

With the unit off, press the "M" button. The display first shows "A", then shows an average of all measurements stored in the unit. **Please note: measurements for each user are averaged and stored separately, be certain that you are viewing the measurements for the correct user.**

Pressing the "M" button again displays the previous value. To view a particular stored memory, press and hold the "M" button to scroll to that stored reading.

Memory full

When the memory has stored 99 results a new measured value is stored by overwriting the oldest value.

Clear all values

If you are sure that you want to permanently remove all stored values, hold down the "M" button (the instrument must have been switched off beforehand) until "CL" appears and then release the button. If you do not want to clear the values, press the On/off button. To permanently clear the memory, press the "M" button while "CL" is flashing.

Individual values cannot be cleared.

5.7. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the START/STOP button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.

5.8. Battery Charge Indicator


Batteries almost discharged

When the batteries are approximately 75% used the battery symbol will flash as soon as the instrument is switched on (if at least one of the batteries still has some charge). Although the instrument will continue to measure reliably, you should obtain replacement batteries.

Batteries discharged— replacements required

When the batteries are discharged, the battery symbol will flash as soon as the instrument is switched on. You cannot take any further measurements and must replace the batteries.

1. Open the battery compartment at the back of the instrument by pushing inwards at the arrow and pulling it out.
2. Replace the batteries – ensure they are correctly connected, as shown on the symbols in the compartment.
3. The memory retains all values although date and time (and possibly also set alarm times) must be reset - the year number therefore flashes automatically after the batteries are replaced.
4. To set date and time, follow the procedure described in Section 4.2.

Which batteries and which procedure?

Use four new, longlife 1.5V AA batteries. Do not use batteries beyond their expiration date. If the monitor is not going to be used for a prolonged period the batteries should be removed.

Using rechargeable batteries

You can also operate this instrument using rechargeable batteries.

- Please use only type "NiMH" reusable batteries!
- If the battery symbol () appears, the batteries must be removed and recharged! They must not remain inside the instrument, as they may become damaged through total discharge even when switched off. The batteries must NOT be discharged in the blood pressure monitor! If you do not intend to use the instrument for a week or more, always remove the rechargeable batteries!
- Recharge these batteries in an external charger and follow manufacturer's instructions carefully.

6. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed. We recommend you switch the device off and take a new measurement making sure you allow time for relaxation. (example: error no. 1)

Error No.

Possible cause(s) / Solutions

ERR 1

The tube may have loosened, or no pulse was detected.

*Ensure cuff connections are tight with proper cuff placement. See section (5.3)

ERR 2

Unnatural pressure impulses. Reason: The arm was moved during the measurement (referred to as an artifact). Repeat measurement keeping arm still.

ERR 3

Repeat measurement keeping arm still. If inflation of the cuff takes too long, the cuff is not correctly seated or the hose connection is not tight. Check connections and repeat.

ERR 5

The difference between systolic and diastolic is excessive. Measure again carefully following proper cuff procedures and ensure measurement under quiet conditions.

HI

The cuff pressure is too high. Relax for 5 minutes and repeat the measurement.*

LO

The pulse is too low (less than 40). Repeat the measurement.*

***If this or any other problem occurs repeatedly, please consult your doctor.**

Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

Malfunction

Remedy

The display remains blank when the instrument is switched on although the batteries are in place.

- Check battery installation.
- If the display is unusual, remove the batteries and then exchange them for new ones. Check polarity.

The pressure does not rise although the pump is running.

- Check the connection of the cuff tube and connect properly.

The device frequently fails to measure, or the values measured are too low or high.

1. Fit the cuff correctly on the arm. (see 5.3)
2. Before starting measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take articles of clothing off if necessary.
3. Measure blood pressure again in complete peace and quiet.

Every measurement results in different values, although the device functions normally and normal values are displayed.

- Please read the following information and the points listed under "Common sources of error". (5.2) Repeat the measurement.
- Blood pressure changes constantly. The observed readings may accurately reflect your pressure.

Blood pressure values differ from those measured by my doctor

- Record daily values and consult your doctor.
- Pressure readings in your doctor's office may be higher due to anxiety.

After the instrument has inflated the cuff the pressure falls very slowly, or not at all. (No reasonable measurement possible).

- Check cuff connections. Ensure the unit has not been tampered with.

7. Care and maintenance

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. **The cuff with bladder must not be washed in a dishwasher, clothes washer, or submerged in water.**
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the monitor!** This invalidates the manufacturer's warranty.
- g) Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

7.1. Accuracy test

We recommend this instrument be tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact Customer Service at 1-800-568-4147 to arrange the test.

8. Warranty

Your Microlife blood pressure monitor is **guaranteed for 5 years** against manufacturer's defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

5 year warranty only applies to the instrument. All accessories including the cuff are guaranteed for one year.

There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by the warranty.

Please note: According to international standards, your monitor should be checked for accuracy every 2 years.

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacture of your blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Essen (RWTüV-Essen).

Never attempt to repair the instrument yourself!

Any unauthorized opening of the instrument invalidates all warranty claims!

9. Certifications

Device standard:

This device is manufactured to meet the European and United States standards for non-invasive blood pressure monitors:

EN1060-1 / 1995

EN1060-3 / 1997

EN1060-4 / 2004

AAMI/ANSI SP10

Electromagnetic compatibility:

Device fulfills the stipulations of the International standard IEC60601-1-2

Clinical testing:

Clinical performance tests were carried out in the U.K. "Association Protocol and British HTP Protocol."

The B.H.S. (British Hypertension Society) clinical protocol was used to measure the accuracy of this product. Blood pressure units using the same measurement technology are graded "AA" for systolic/diastolic accuracy by independent investigators using the BHS protocol. This is the highest grading available for blood pressure monitors. Please see bhsoc.org for more information. (Uses the same algorithm as B.H.S. graded "AA" model number 3BT0-1)

US patents:

6,485,429 (Irregular Heartbeat Feature)

6,447,457 (Averaging Mode)

10. Technical specifications:

Weight:	560 g (with cuff and batteries)
Size:	149 (W) x 96 (L) x 39 (H) mm
Storage temperature:	-20 to +50°C
Humidity:	15 to 90% relative humidity maximum
Operation temperature:	10 to 40°C (50°F - 104°F)
Display:	LCD-Display (Liquid Crystal Display)
Measuring method:	oscillometric
Pressure sensor:	capacitive
Measuring range:	
SYS/DIA:	30 to 280 mmHg
Pulse:	40 to 200 per minute
Cuff pressure display range:	0-299 mmHg
Memory:	Automatically stores up to 198 measurement readings and 99 sets of memory
Measuring resolution:	1 mmHg
Accuracy:	pressure within ± 3 mmHg pulse ± 5 % of the reading
Power source:	a) 4 AA batteries, 1.5 V b) AC adapter 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC, use only AC Adapter Model # HMA-AC2)
Accessories:	cuff type S101 for arm circumference 24 - 33.7 cm (9.5" - 13.25") included cuff type S102 for arm circumference 30.5-40.64 cm (12" - 16") sold as optional accessory

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Microlife USA, Inc.