# PREMIUM AUTOMATIC Blood Pressure Monitor

Product #344532  
Instruction Manual  
Instrucciones de Uso

Questions? Preguntas?  
1-866-464-6184

PLEASE NOTE:  
THIS MEDICAL INSTRUMENT MUST BE USED ACCORDING  
TO INSTRUCTIONS TO ENSURE ACCURATE READINGS.

# Automatic Blood Pressure Monitor - Instruction Manual

## Table of contents

---

### **1. Introduction**

- 1.1. Features of your blood pressure monitor
- 1.2. Important information about self measurement

### **2. Important information on blood pressure and its measurement**

---

- 2.1. How does high/low blood pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done if regular high/low values are obtained?

### **3. Components of your blood pressure monitor**

---

### **4. Using your monitor for the first time**

---

- 4.1. Activating the pre-installed batteries
- 4.2. Setting the date and time
- 4.3. Cuff tube connection
- 4.4. Select the user

### **5. Measurement Procedure**

---

- 5.1. Before measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Irregular Heartbeat Detector
- 5.6. Memory - Displaying the last 99 measurements
- 5.7. Discontinuing a measurement
- 5.8. Battery Charge Indicator
- 5.9. Using an AC Adapter

### **6. Error Messages/Troubleshooting**

---

### **7. Care and maintenance**

---

- 7.1. Accuracy Test

### **8. Warranty**

---

### **9. Certifications**

---

### **10. Technical specifications**

---

### **11. How to contact us**

## 1. Introduction

---

### 1.1. Features of your CVS automatic blood pressure monitor

Your CVS blood pressure monitor with Irregular Heartbeat Detection is a fully automatic digital blood pressure measuring device for use by adults on the upper arm at home or in your doctor's/nurse's office. It enables very fast and reliable measurement of systolic and diastolic blood pressure as well as pulse by way of the oscillometric method. This device offers clinically proven accuracy and has been designed to be user friendly.

Before using, please read this instruction manual carefully and then keep it in a safe place. Please contact your doctor for further questions on the subject of blood pressure and its measurement.

### Attention!

#### 1.2. Important information about self-measurement

- **Self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.**
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In cases of irregular heartbeat (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with your doctor.

#### Electromagnetic interference:

The device contains sensitive electronic components. Avoid strong electrical or electromagnetic fields in the direct vicinity of the device (e.g. mobile telephones, microwave ovens). These can lead to temporary impairment of the measuring accuracy.

## 2. Important information on blood pressure and its measurement

---

### 2.1. How does high/low blood pressure arise?

Your level of blood pressure is determined in the circulatory center of the brain and adjusts to a variety of situations through feedback from the nervous system. To adjust blood pressure, the strength and frequency of the heart (Pulse), as well as the width of circulatory blood vessels is altered. Blood vessel width is effected by fine muscles in the blood vessel walls.

Your level of arterial blood pressure changes periodically during heart activity: During the “blood ejection” (Systole) the value is highest (systolic blood pressure value). At the end of the heart’s “rest period” (Diastole) pressure is lowest (diastolic blood pressure value).

Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

### 2.2. Which values are normal?

Blood pressure is very high if your diastolic pressure is above 90 mmHg and/or your systolic blood pressure is over 160 mmHg, **while at rest**. In this case, please consult your physician immediately. Long-term values at this level endanger your health due to continual damage to the blood vessels in your body.

If your systolic blood pressure values are between 140 mmHg and 159mmHg and/or the diastolic blood pressure values between 90 mmHg and 99mmHg, consult your physician. Regular self-checks are necessary.

If you have blood pressure values that are too low, (i.e., systolic values under 105mmHg and/or diastolic values under 60 mmHg), consult your physician.

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. You can detect possible changes in your values early and react appropriately.

If you are undergoing medical treatment to control your blood pressure, keep a record of values along with time of day and date. Show these values to your physician. **Never use the results of your measurements to independently alter the drug doses prescribed by your physician.**

#### Which values are normal?

The following standards for assessing high blood pressure (without regard to age) have been established by the National Institutes of Health JNCV7, 2003.

Category	Systolic (mmHg)	Diastolic (mmHg)
Normal	<120	<80
Pre-Hypertension	120-139	80-89
<b>Hypertension</b>		
Stage 1 Hypertension	140-159	90-99
Stage 2 Hypertension	>160	>100

### **Further information**

- If your values are mostly normal under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called “labile hypertension.” Consult your doctor.
- **Correctly measured diastolic blood pressure values above 120mmHg and require immediate medical treatment.**

### **2.3. What can be done if regular high or low values are obtained?**

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organ malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take to reduce and even prevent high blood pressure. These measures must be permanent lifestyle changes.

#### **A) Eating habits**

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

#### **B) Previous illnesses**

- Consistently follow all medical instructions for treating illness such as:
  - Diabetes (Diabetes mellitus or sugar diabetes)
  - Fat metabolism disorder
  - Gout

#### **C) Habits**

- Give up smoking completely.
- Drink only moderate amounts of alcohol.
- Restrict your caffeine consumption (coffee, tea, chocolate).


## D) Your Physical condition

- **After a medical examination, and with your doctor's approval and direction; exercise.**
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your exercise routine. You must receive advise regarding the type and extent of exercise that is appropriate for you.

## 3. Components of your blood pressure monitor

---

### a) Measuring unit


### Replacing the removable card

You can replace the slot-in card by pulling it out to the side and replacing the paper insert. It may be helpful to have your doctor note down the medication dosage or an emergency telephone number on the card.

## b) Upper arm cuffs:

Type S101 for arm circumference 24–33.7 cm or  
9.5" - 13.25"

Type S102 for arm circumference 30.5 - 40.64 cm  
12" - 16" (available for purchase  
separately. CVS model # 271237.  
Use only CVS or Microlife brand cuffs!)

Please Note:

Arm circumference should be measured with a measuring tape in the middle of the relaxed upper arm. Do not force cuff connection into the opening. Make sure the cuff connection is not pushed into the AC adapter port. If the cuff is too small, call 1-866-464-6184 for further information. You may prefer to use one of our wrist type monitors.

## 4. Using your monitor for the first time

---

### 4.1 Activating the pre-installed batteries

Pull out the protective strips projecting from the battery compartment.

### 4.2 Setting the date and time

1. After the new batteries are activated, the year number flashes in the display. You can advance the year by pressing the "M" button. To confirm and then set the month, press the time button.
2. You can now set the month using the «M» button. To confirm and then set the day, press the time button.
3. Please follow the instructions above to set the day, hour and minutes.
4. Once you have set the last minute and pressed the time button, the date and time are set and the time is displayed.
5. If you want to change the date and time, press and hold the time button down for approx. 3 seconds until the year number starts to flash. Now you can enter the new values as described above.

### **4.3. Cuff tube connection**

Insert the cuff tube into the opening on the side of the monitor.

### **4.4. Select the User**

This blood pressure monitor is designed to store 99 measurements for each of two users. Before taking a measurement, be certain that the correct user has been selected.

- a) With the unit off, press and hold the TIME button for 3 seconds until the user icon in the upper left corner of the LCD screen flashes.
- b) Press the "M" button to toggle between users.
- c) Press the START button to make your selection.

## 5. Measurement Procedure

---

Please note: You should always be seated before and during measurement.

### 5.1. Before measurement:

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always compare measurements taken at the same time of day, since blood pressure changes during the course of the day.

### 5.2. Common sources of error:

**Note: Comparable blood pressure measurements always require the same conditions! Conditions should always be quiet.**

- All efforts by the user to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not flex any of the muscles in the measurement arm during the measurement. Use a cushion for support if necessary.
- If the arm artery lies considerably lower or higher than the heart, an erroneously high or low blood pressure will be measured! Each 15cm (6 inch) difference in height between your heart and the cuff results in a measurement error of 10 mmHg!
- Cuffs that are too narrow or too short result in false measurement values. Selecting the correct cuff is extremely important. Cuff size is dependent upon the circumference of the arm (measured in the center). The permissible range is printed on the cuff. If this is not suitable for your use, please call 1-866-464-6184.

**Note: Only use clinically approved Cuffs!**


- A loose cuff or a sideways protruding air pocket causes false measurement values.
- With repeated measurements, blood accumulates in the arm, which can lead to false results. Consecutive blood pressure measurements should be repeated after a 1 minute pause or after your arm has been held up in order to allow the accumulated blood to flow away.

### 5.3. Fitting the cuff


- Pass the end of the cuff through the flat metal ring so that a loop is formed. The velcro closer must be facing outward. (Ignore this step if the cuff has already been prepared.)
- Place the cuff over the left upper arm so that the tube is closer to your lower arm.
- Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately  $3/4''$  to  $1''$  (2 to 3cm) above the elbow and that the tube is closer to the inner side of the arm.

**Important!** The **small white arrow (Artery Mark)** on the cuff must lie exactly over the artery which runs down the inner side of the arm.

- Tighten the cuff by pulling the end and close the cuff by affixing the velcro.
- There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed. Cuffs that don't fit properly result in false measurement values. Measure your arm circumference if you are not sure of proper fit. Larger cuffs are available (Section 3B).
- Lay your arm on a table (palm upward) so the cuff is at the same height as your heart. Make sure the tube is not kinked.
- Remain seated quietly for at least two minutes before you begin the measurement.

**ARTERY MARK**

## Comment:

If it is not possible to fit the cuff to your left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.

Comparable blood pressure measurements always require the same conditions (Relax for several minutes before a reading).

Cuff on  
right arm

## 5.4. Measuring procedure

After the cuff has been appropriately positioned the measurement can begin:

- Press the ON/OFF/START button. The pump begins to inflate the cuff. In the display, the increasing cuff pressure is continually displayed.
- After automatically reaching an individual pressure, the pump stops and the pressure slowly falls. The cuff pressure is displayed during the measurement.
- When the device has detected your pulse, the heart symbol in the display begins to blink.
- When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values, as well as the pulse are now displayed and the blue backlight is illuminated.
- The appearance of this symbol  signifies that an irregular heartbeat was detected. This indicator is only a caution. It is important that you be relaxed, remain still and do not talk during measurements. NOTE: We recommend contacting your physician if you see this indicator frequently.
- The measurement results are displayed until you switch the device off. If no button is pressed for 1 minute, the device switches off automatically.

Pumping  
Pressure

A digital display showing the number 162 in a large, black, seven-segment font.

Measuring

A digital display showing the number 146 in a large, black, seven-segment font. A small heart symbol is positioned below the number.

Systolic

A digital display showing the number 129 in a large, black, seven-segment font.

Diastolic

A digital display showing the number 78 in a large, black, seven-segment font.

Pulse

A digital display showing the number 5 in a large, black, seven-segment font.

Measurement  
complete

## 5.5. Irregular Heartbeat Detector


This symbol  indicates that certain pulse irregularities were detected during the measurement. In this case, the result may deviate from your normal basal blood pressure – repeat the measurement. In most cases, this is no cause for concern. However, if the symbol appears on a regular basis (e.g. several times a week with measurements taken daily), we advise you to tell your doctor. Please show your doctor the following explanation:

### Information for the doctor on frequent appearance of the Irregular Heartbeat Symbol

This instrument is an oscillometric blood pressure monitor device that also analyses pulse frequency during measurement. The instrument is clinically tested.

If pulse irregularities occur during measurement, the irregular heartbeat symbol is displayed after the measurement. If the symbol appears more frequently (e.g. several times per week on measurements performed daily) or if it suddenly appears more often than usual, we recommend the patient to seek medical advice.

The instrument does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

## 5.6. Memory – displaying the last 99 measurements

At the end of a measurement, this monitor automatically stores each result with date and time. This unit stores 99 memories for each of 2 users.

### Viewing the stored values

With the unit off, press the "M" button. The display first shows "A", then shows an average of all measurements stored in the unit. **Please note: measurements for each user are averaged and stored separately, be certain that you are viewing the measurements for the correct user.**

Pressing the "M" button again displays the previous value. To view a particular stored memory, press and hold the "M" button to scroll to that stored reading.

### Memory full

When the memory has stored 99 results a new measured value is stored by overwriting the oldest value.

### Clear all values

If you are sure that you want to permanently remove all stored values, hold down the "M" button (the instrument must have been switched off beforehand) until "CL" appears and then release the button. If you do not want to clear the values, press the On/off button. To permanently clear the memory, press the "M" button while "CL" is flashing.

Individual values cannot be cleared.

## 5.7. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason (e.g the patient feels unwell), the ON/OFF button can be pressed at any time. The device then immediately lowers the cuff pressure automatically.


**Start**  
ON/OFF

## 5.8. Battery Charge Indicator

---


### Batteries almost discharged

When the batteries are approximately 75% used the battery symbol will flash as soon as the instrument is switched on (if at least one of the batteries still has some charge). Although the instrument will continue to measure reliably, you should obtain replacement batteries.

### Batteries discharged— replacements required

When the batteries are discharged, the battery symbol will flash as soon as the instrument is switched on. You cannot take any further measurements and must replace the batteries.

1. Open the battery compartment at the back of the instrument by pushing inwards at the two arrows and pulling it out.
2. Replace the batteries – ensure they are correctly connected, as shown on the symbols in the compartment.
3. The memory retains all values although date and time (and possibly also set alarm times) must be reset - the year number therefore flashes automatically after the batteries are replaced.
4. To set date and time, follow the procedure described in Section 4.2.

### Which batteries and which procedure?

Use four new, longlife 1.5V AA batteries. Do not use batteries beyond their expiration date. If the monitor is not going to be used for a prolonged period the batteries should be removed.

### Using rechargeable batteries

You can also operate this instrument using rechargeable batteries.

- Please use only type "NIMH" reusable batteries!
- If the battery symbol (  ) appears, the batteries must be removed and recharged! They must not remain inside the instrument, as they may become damaged through total discharge even when switched off. The batteries must NOT be discharged in the blood pressure monitor! If you do not intend to use the instrument for a week or more, always remove the rechargeable batteries!
- Recharge these batteries in an external charger and follow manufacturer's instructions carefully.

## 5.9. Using an AC Adapter

---

You may also operate this monitor using an AC Adapter.

Call 1-866-464-6184 for more information or to purchase an AC Adapter.


1. Ensure that the AC adapter and cable are not damaged.
2. Plug the adapter cable into the AC Adapter port on the blood pressure monitor.
3. Plug the adapter into your electrical outlet.

When the AC adapter is connected, no battery current is consumed.

**AC Adapter sold separately. Use only CVS brand AC Adapter (Model #306409) or Microlife brand AC Adapters.**

## 6. Error Messages/Troubleshooting

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed. We recommend you switch the device off and take a new measurement making sure you allow time for relaxation. (example: error no. 1)

<b>Error No.</b>	<b>Possible cause(s) / Solutions</b>
ERR 1	The tube may have loosened, or no pulse was detected. *Ensure cuff connections are tight with proper cuff placement. See section (5.3)
ERR 2	Unnatural pressure impulses. Reason: The arm was moved during the measurement (Artefact). Repeat measurement keeping arm still.
ERR 3	Repeat measurement keeping arm still. If inflation of the cuff takes too long, the cuff is not correctly seated or the hose connection's not tight. Check connections and repeat.
ERR 5	The difference between systolic and diastolic is excessive. Measure again carefully following proper cuff procedures and ensure measurement under quiet conditions.
HI	The cuff pressure is too high. Relax for 5 minutes and repeat the measurement.*
LO	The pulse is too low (less than 40). Repeat the measurement.*

**\*If this or any other problem occurs repeatedly, please consult your doctor.**

## Other possible errors and their solutions

If problems occur when using the device, the following points should be checked:

<b>Malfunction</b>	<b>Remedy</b>
The display remains blank when the instrument is switched on although the batteries are in place.	<ul style="list-style-type: none"><li>• Check battery installation.</li><li>• If the display is unusual, remove the batteries and then exchange them for new ones. Check polarity.</li></ul>
The pressure does not rise although the pump is running.	<ul style="list-style-type: none"><li>• Check the connection of the cuff tube and connect properly.</li></ul>
The device frequently fails to measure, or the values measured are too low or high.	<ol style="list-style-type: none"><li>1. Fit the cuff correctly on the arm. (see 5.3)</li><li>2. Before starting measurement make sure that the cuff is not too tight and that clothing is not exerting pressure on the arm. Take articles of clothing off if necessary.</li><li>3. Measure blood pressure again in complete peace and quiet.</li></ol>
Every measurement results in different values, although the device functions normally and normal values are displayed.	<ul style="list-style-type: none"><li>• Please read the following information and the points listed under "Common sources of error". (5.2) Repeat the measurement.</li><li>• Blood pressure changes constantly. The observed readings may accurately reflect your pressure.</li></ul>
Blood pressure values differ from those measured by my doctor	<ul style="list-style-type: none"><li>• Record daily values and consult your doctor.</li><li>• Pressure readings in your doctor's office may be higher due to anxiety.</li></ul>
After the instrument has inflated the cuff the pressure falls very slowly, or not at all. (No reasonable measurement possible).	<ul style="list-style-type: none"><li>• Check cuff connections. Ensure the unit has not been tampered with.</li></ul>

## 7. Care and maintenance

---

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this cuff carefully and avoid all types of stress through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. **The cuff with bladder (see 8.1 below) must not be washed in a dishwasher, clothes washer, or submerged in water.**
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the monitor!** This invalidates the manufacturer's warranty.
- g) Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

### 7.1. Accuracy test

We recommend this instrument is tested for accuracy every two years or after mechanical impact (e.g. being dropped). Please contact Customer Service at 1-866-464-6184 to arrange the test.

## 8. Warranty

---

Your CVS blood pressure monitor is **guaranteed for a lifetime** against manufacturers' defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

Lifetime warranty only applies to the instrument. All accessories including the cuff are guaranteed for one year.

There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by the warranty.

Please note: According to international standards, your monitor should be checked for calibration every 2 years.

In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacture of your blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Essen (RWTüV-Essen).

**Never attempt to repair the instrument yourself!**

Any unauthorized opening of the instrument invalidates all warranty claims!

## 9. Certifications

---

### **Device standard:**

This device is manufactured to meet the European and United States standards for non-invasive blood pressure monitors:  
EN1060-1 / 1995  
EN1060-3 / 1997

### **Electromagnetic compatibility:**

DIN 58130, NIBP – clinical investigation  
ANSI / AAMI SP10, NIBP – requirements

### **Clinical testing:**

Device fulfills the stipulations of the International standard EN 60601-1-2  
Clinical performance tests were carried out in the U.K. “Association Protocol and British HTP Protocol.”

The B.H.S. (British Hypertension Society) clinical protocol was used to measure the accuracy of this product. Blood pressure units using the same measurement technology are graded “AA” for systolic/diastolic accuracy by independent investigators using the BHS protocol. This is the highest grading available for blood pressure monitors. Please see [bhsoc.org](http://bhsoc.org) for more information. (Uses the same algorithm as B.H.S. graded “AA” model number 3BT0-1)

US patents:  
6,485,429 (Irregular Heartbeat Feature)

## 10. Technical specifications:

<b>Weight:</b>	505 g (with batteries)
<b>Size:</b>	125 (W) x 151 (L) x 95 (H) mm
<b>Storage temperature:</b>	-5 to +50°C (-4°F - 122°F)
<b>Humidity:</b>	15 to 90% relative humidity maximum
<b>Operation temperature:</b>	10 to 40°C ( 50°F - 104°F)
<b>Display:</b>	LCD-Display (Liquid Crystal Display)
<b>Measuring method:</b>	oscillometric
<b>Pressure sensor:</b>	capacitive
<b>Measuring range:</b>	
<b>SYS/DIA:</b>	30 to 280 mmHg
<b>Pulse:</b>	40 to 200 per minute
<b>Cuff pressure display range:</b>	0–299 mmHg
<b>Memory:</b>	Automatically stores the last 99 measurements for 2 users (total 198)
<b>Measuring resolution:</b>	1 mmHg
<b>Accuracy:</b>	pressure within $\pm 3$ mmHg pulse $\pm 5$ % of the reading
<b>Power source:</b>	a) 4 AA batteries, 1.5 V b) AC adapter 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC)
<b>Accessories:</b>	cuff type S101 for arm circumference 24 - 33.7 cm (9.5" - 13.25") included cuff type S102 for arm circumference 30.5–40.64 cm (12" - 16") sold as optional accessory

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**©2004 CVS/pharmacy**
CVS Pharmacy, Inc.