# MANUAL INFLATE

## Blood Pressure Monitor

Product #271243  
Instruction Manual  
Instrucciones de Uso

Questions? Preguntas?  
1-866-464-6184

PLEASE NOTE:  
THIS MEDICAL INSTRUMENT MUST BE USED ACCORDING  
TO INSTRUCTIONS TO ENSURE ACCURATE READINGS.

# CVS Semi-automatic Blood Pressure Monitor

## Instruction Manual

### **Table of contents**

---

#### **1. Introduction**

---

- 1.1 Features of your blood pressure monitor.
- 1.2 Important information about self-measurement

#### **2. Important information on the subject of blood pressure and its measurement**

---

- 2.1. How does high/low blood-pressure arise?
- 2.2. Which values are normal?
- 2.3. What can be done, if regular high/low values are obtained?

#### **3. Components of your blood pressure monitor**

---

#### **4. Putting your blood pressure monitor into operation**

---

- 4.1. Inserting the batteries
- 4.2. Tube connection
- 4.3. Setting the Time and Date

#### **5. Measurement procedure**

---

- 5.1. Before a measurement
- 5.2. Common sources of error
- 5.3. Fitting the cuff
- 5.4. Measuring procedure
- 5.5. Inflating further
- 5.6. Discontinuing a measurement
- 5.7. Reading recall - Displaying the last measurement

#### **6. Error messages/troubleshooting**

---

#### **7. Care and maintenance**

---

#### **8. Warranty**

---

#### **9. Certifications**

---

#### **10. Technical specifications**

---

#### **11. How to contact us**

---

## 1. Introduction

---

### 1.1. Features of your CVS semi-automatic blood pressure monitor, model #271243

Your blood pressure monitor is a semi-automatic, digital blood pressure measuring device for use on the upper arm. It enables very fast and reliable measurement of the systolic and diastolic blood pressure as well as the pulse through the oscillometric method of measuring.

The device offers a clinically tested measurement accuracy and has been designed to be user-friendly.

The large display allows both operational status and cuff pressure to be visible for the entire duration of the measurement.

Before using, please read this instruction manual carefully and then keep it in a safe place.

For further questions on the subject of blood pressure and its measurement, please contact your doctor.

#### **Attention!**

### 1.2. Important information about self-measurement

- **Self-measurement means Control**, not diagnosis or treatment. Unusual values must always be discussed with your doctor. **Under no circumstances should you alter the dosages of any drugs prescribed by your doctor.**
- The pulse display is **not** suitable for checking the frequency of heart pacemakers!
- In cases of irregular heartbeat (Arrhythmia), measurements made with this instrument should only be evaluated after consultation with your doctor.

#### **Electromagnetic interference:**

The device contains sensitive electronic components (Microcomputer). Therefore, avoid strong electrical or electromagnetic fields in the direct vicinity of the device, such as cellular telephones and microwave ovens. These devices, when in operation close to your blood pressure monitor, can lead to temporary impairment of the measuring accuracy.

## 2. Important information on blood pressure and its measurement

---

### 2.1. How does high/low blood pressure arise?

Your level of blood pressure is determined in a part of the brain called the circulatory center. This part of your brain responds to a situation through feedback with the nervous system. To adjust the blood pressure, the strength and frequency of the heart (Pulse) as well as the width of circulatory blood vessels is altered. Blood vessel width is altered by fine muscles in the blood vessel walls.

The level of arterial blood pressure changes periodically during heart activity: During blood ejection from the heart (Systole) the value is maximal (systolic blood pressure value). At the end of the heart's 'rest period', Diastole or diastolic blood pressure value, blood pressure falls to a low value. Blood pressure values must lie within certain normal ranges in order to prevent particular diseases.

## 2.2 Which values are normal?

Blood pressure is too high if at rest, your diastolic pressure is above 90mmHg and/or the systolic blood pressure is over 160mmHg. **You should bring these values to the attention of your doctor immediately. Long-term values at this level endanger your health due to the associated advancing damage to the blood vessels in your body.**

**Should the systolic blood pressure values lie between 140mmHg and 160mmHg and/or the diastolic blood pressure values between 90mmHg and 95mmHg, please consult your doctor, and check your blood pressure regularly.**

**With blood pressure values that are too low, i.e. systolic values under 105mmHg and/or diastolic values under 60mmHg, please consult your doctor.**

Even with normal blood pressure values, a regular self-check with your blood pressure monitor is recommended. In this way you can detect possible changes in your values early and react appropriately, in consultation with your physician.

If you are undergoing medical treatment to control your blood pressure, please keep a record of your blood pressure levels by carrying out regular self-measurements at specific times of the day. Show these values to your doctor. **Never use the results of your measurements to discontinue or independently alter the drug doses prescribed by your doctor.**

The following standards for assessing high blood pressure (without regard to age) have been established by the National Institutes of Health JNCVI.

<b>Category</b>	<b>Systolic (mmHg)</b>	<b>Diastolic (mmHg)</b>
Optimal	<120	<80
Normal	<130	<85
High Normal	130-139	85-89
<b>Hypertension</b>		
Stage 1	140-159	90-99
Stage 2	160-179	100-109
Stage 3	≥180	≥110

### Further information

- If your values are mostly standard under resting conditions but exceptionally high under conditions of physical or psychological stress, it is possible that you are suffering from so-called "labile hypertension". Please consult your doctor if you suspect that this might be the case.

### **Warning:**

- **Diastolic blood pressure values above 120mmHg require immediate medical treatment.**

### 2.3. What can be done if regular high or low values are obtained?

- a) Consult your doctor.
- b) Increased blood pressure values (various forms of hypertension) are associated with considerable health risks over time. Arterial blood vessels in your body are endangered due to constriction caused by deposits in the vessel walls (Arteriosclerosis). A deficient supply of blood to important organs (heart, brain, muscles) can result from arteriosclerosis. Furthermore, the heart will become structurally damaged with increased blood pressure values.
- c) There are many different causes of high blood pressure. We differentiate between the common primary (essential) hypertension, and secondary hypertension. The latter group can be ascribed to specific organic malfunctions. Please consult your doctor for information about the possible origins of your own increased blood pressure values.
- d) There are measures which you can take to reduce and prevent high blood pressure. These measures must be permanent lifestyle changes.

#### A) Eating habits

- Strive for a normal weight corresponding to your age. See your doctor for your ideal weight.
- Avoid excessive consumption of common salt.
- Avoid fatty foods.

#### B) Previous illnesses

- Consistently follow all medical instructions for treating illnesses such as:
- Diabetes (Diabetes mellitus or sugar diabetes)
- Fat metabolism disorder
- Gout

#### C) Habits


- Give up smoking completely
- Drink only moderate amounts of alcohol
- Restrict your caffeine consumption (coffee, tea, chocolate)

#### D) Physical constitution

- After a preliminary medical examination, exercise regularly.
- Choose sports which require stamina and avoid those which require strength.
- Avoid reaching the limit of your performance.
- With previous illnesses and/or an age of over 40 years, please consult your doctor before beginning your sporting activities. He/She will advise you regarding the type and extent of types of sport that are possible for you.

### 3. Components of the blood-pressure monitor

#### a) Measuring unit

#### b) Upper arm cuffs:

The cuff is for arm circumference

24–34 cm or 9.5" - 13.25" (included)

CVS sells a large cuff for arm circumference 31-41 cm or

12" - 16" (sold separately in all CVS stores)

### 4. Operation of your blood pressure monitor.


#### 4.1. Inserting the batteries

After you have unpacked your device, insert the batteries. The battery compartment is located on the underside of the device (see illustration).

a) Remove the battery cover.

b) Insert the batteries (4 x size AA 1.5 V), following the indicated polarity.

c) If the battery warning appears in the display, the batteries are discharged and must be replaced with new ones.

#### Attention!

- After the battery warning appears, the device is blocked until the batteries have been replaced.
- Please use «AA» Long-Life or Alkaline 1.5 V Batteries. The

Low battery indicator


use of rechargeable batteries is not recommended.

- If your blood pressure monitor is not used for long periods, please remove the batteries from the device.

### **Functional check:**

Hold on the on/off button down to test all the display elements.

When functioning correctly icons will appear.

## **4.2. Tube connection**

### **• Cuff**

Set up (or check) the tube connections between cuff, instrument and pump ball as shown in the diagram.

## **4.3. Setting the Time and Date**

To set the time and date press the Set Time/Date button and hold it down for five seconds. When you let go the number you'll see a flashing four digit number (probably 2001 or 2002 to 2006). This number represents the year. Use the MEMORY key to scroll up (you cannot scroll backwards so scroll slowly) to the correct year. When you reach the correct year stop scrolling and press the Set Time/Date button to set the year.

After doing this the number flashing will be the month. Again, use the MEMORY key to scroll to the appropriate month and stop scrolling then press the Set Time/Date button to set the month and move on to setting the day. Use the MEMORY key to scroll up to find the correct day of the month and press the Set Time/Date button to set it.

Now that you've set the correct date, the number flashing will represent the time. The first two digits are the hour. Use the MEMORY button to scroll through the hours. Note that for morning hours before noon, there is no AM displayed but after noon PM will display. Again, scroll through the time using the MEMORY button to scroll and the Set Time/Date button to set the time.

## **5. Measurement Procedure**

---

Please note: You should always be seated before and during measurement

### **5.1. Before the measurement:**

- Avoid eating and smoking as well as all forms of exertion directly before measurement. These factors influence the measurement result. Find time to relax by sitting in an armchair in a quiet atmosphere for about ten minutes before measurement.
- Remove any garment that fits closely to your upper arm.
- Always measure on the same arm (normally left).
- Always perform measurements at the same time of day, since the blood pressure changes during the course of the day.

## 5.2. Common sources of error:

### Note:

**Comparable blood pressure measurements always require the same conditions!**

**These are normally quiet conditions.**

- Efforts by the patient to support the arm can increase blood pressure. Make sure you are in a comfortable, relaxed position and do not activate any of the muscles in the measurement arm during the measurement. Use a cushion for support and rest your arm on a table.
- If the upper arm artery lies considerably lower or higher than the heart, a false high or low blood pressure will be measured!

**\*A variation of 6" between the cuff and heart level can result in a reading error of + or - 10mmHg.**

- Cuffs that do not fit properly result in false measurement values. Selecting the correct cuff is extremely important. The cuff size is dependent upon the circumference or distance around your upper arm measured in the center. The permissible range is printed on the cuff. If this is not suitable for your use, please contact CVS.

**Note:** Only use clinically approved **Microlife or CVS Cuff Model # 271237.**


- A loose cuff or a sideways protruding air-pocket causes false measurement values.
- Measurements should be done after a 5 minute rest to ensure accuracy.
- Consecutive blood pressure measurements should be performed after a 1 minute pause.

## 5.3. Fitting the cuff


- a) Pass the end of the cuff through the flat metal ring so that a loop is formed. The Velcro closer must be facing outwards. (Ignore this step if the cuff has already been prepared.)


- b) Place the cuff over the left upper arm so that the tube is closer to your lower arm. Align the cuff artery mark (arrow) over the arm artery than runs down the inner side of the arm.

- c) Lay the cuff on the arm as illustrated. Make certain that the lower edge of the cuff lies approximately 3/4" to 1" (2 to 3 cm) above the elbow and that the tube is closer to the inner side of the arm.

- d) Tighten the cuff by pulling the end and close the cuff by affixing the velcro. Make sure that the cuff Index end falls between the OK lines when the cuff is in po-


sition. If it doesn't, you are using the wrong cuff size.

- e) There should be little free space between the arm and the cuff. You should be able to fit 2 fingers between your arm and the cuff. Clothing must not restrict the arm. Any piece of clothing which does must be removed.
- f) Lay your arm on a table (palm upward) so the cuff is at the same height as the heart. Make sure that the tube is not kinked.
- g) Remain seated quietly for five minutes before you begin the measurement.

**Comment:**


If it is not possible to fit the cuff to the left arm, it can also be placed on your right arm. However, all measurements should be made using the same arm.


**5.4. Measuring procedure**

After the cuff has been appropriately positioned, the measurement can begin:


- a) Switch on the instrument by pressing the On/Off/START button. Many icons appear in the display for 2 seconds: (Diag.), A short triple beep follows to indicate that the instrument is ready. A «0» is shown in the display.
- b) Grasp the pump ball with your free hand and Inflate the cuff to a pressure at least 40 mmHg higher than the expected systolic pressure. If you do not know the value of this pressure, inflate to a value of approximately 160–180 mmHg. The cuff pressure is shown continuously on the display.

- c) After the inflation process lay down the pump ball and remain in a quiet and relaxed sitting position. The measurement now proceeds on its own.  
If no measurement takes place and all that appears in the display is a flashing, upwards pointing arrow, pump the ball once again and inflate to a higher pressure «Inflate further». When the instrument detects a pulse, a heart symbol in the display starts to flash and a beep is heard for every heart beat.


- d) When the measurement has been concluded, a long beep tone sounds. The measured systolic and diastolic blood pressure values as well as the pulse are now displayed. The cuff pressure can now be quickly released by pressing the release valve at the pump ball.

### 5.5. Inflating further


If the cuff is not inflated sufficiently, the measurement is stopped after a few seconds and a flashing arrow pointing upwards appears in the display. It is then necessary to inflate at least 20 mmHg higher than the previous value. The instruction to inflate further can appear several times if the pressure is still not sufficient.

### Attention:


If the cuff is inflated too far (over 300 mmHg), «HI» appears in the display and a warning signal is sounded. The cuff pressure must then be immediately released using the quick-release valve at the pump ball!

As soon as the cuff is inflated sufficiently and the measurement is proceeding, continue as described under points c + d.


### 5.6. Discontinuing a measurement

If it is necessary to interrupt a blood pressure measurement for any reason, press the quick release valve at the pump ball. The device immediately lowers the cuff-pressure automatically.

### 5.7. Reading Recall - Displaying the last measurement

The measured results are stored in the instrument. This unit stores 30 measurements. You can call up the stored values at any time by pressing the MEMORY button. The first memory displayed will be the most recent measurement taken. Keep pushing the button to scroll back through memories.

#### Further information


Measurements should not occur soon after each other, to ensure accurate readings. Wait at least **1 minute** in a relaxed position, sitting or lying, before you repeat a measurement.

## 6. Error messages/Troubleshooting

---

If an error occurs during a measurement, the measurement is discontinued and a corresponding error code is displayed.

(example: error no. 1)


Error No.	Possible cause(s)
ER 1	The systolic pressure was determined but afterwards the cuff pressure fell below 20 mmHg. The tube may have become unplugged. Further possible cause: No pulse was detected.
ER 2	The pressure dropped below 20 mmHg without a measurement taking place. The Release valve was pressed or the arm moved during the measurement.
ER 5	Unit detected an unacceptable difference between systolic and diastolic pressures. Take another reading. Contact your doctor if this continues.

## Other possible malfunctions and possible solutions:

<b>Malfunction</b>	<b>Remedy</b>
The display remains blank when the instrument is switched on although the batteries are in place.	<ol style="list-style-type: none"><li>1. Check battery polarity and correct if necessary.</li><li>2. If the display is unusual, remove the batteries and insert new batteries.</li></ol>
The pressure does not rise although you are pumping the inflation ball.	<ul style="list-style-type: none"><li>• Check the connection of the cuff tube and connect properly if necessary.</li></ul>
The device frequently fails to measure the blood pressure values, or the values measured are too low (too high).	<ol style="list-style-type: none"><li>1. Fit the cuff correctly on the arm.</li><li>2. Before starting a measurement make sure that the cuff is not fitted too tightly. Ensure clothing is not exerting pressure on the arm above the measuring position. Take off articles of clothing if necessary.</li><li>3. Measure blood pressure again in complete peace and quiet.</li></ol>
Every measurement results in different values, although the device functions normally and normal values are displayed.	<ul style="list-style-type: none"><li>• Blood pressure changes continually. Your readings may be accurate. Consult your doctor.</li></ul>
The blood pressure values measured differ from those measured by the doctor	<ul style="list-style-type: none"><li>• Record the daily development of the values and consult your doctor about them.</li><li>• Anxiety during a doctor visit can result in higher readings.</li></ul>

## 7. Care and maintenance

- a) Do not expose the device to either extreme temperatures, humidity, dust or direct sunlight.
- b) The cuff contains a sensitive air-tight bubble. Handle this carefully and avoid all types of straining through twisting or buckling.
- c) Clean the device with a soft, dry cloth. Do not use gas, thinners or similar solvents. Spots on the cuff can be removed carefully with a damp cloth and soapsuds. **The cuff must not be washed in a dishwasher, clothes washer, or submerged in water.**
- d) Handle the tube carefully. Do not pull on it. Do not allow the tubing to kink and keep it away from sharp edges.
- e) Do not drop the monitor or treat it roughly in any way. Avoid strong vibrations.
- f) **Never open the monitor!** This invalidates the manufacturer's warranty.

### Further information

Blood pressure is subject to fluctuations even in healthy people.

**Comparable measurements always require the same conditions (quiet conditions)! If fluctuations in readings are larger than 15mmHg, and/or you hear irregular pulse tones, consult your doctor.** In order to receive market clearance from governmental bodies, this device was subjected to strict clinical tests. The computer program used to measure blood pressure values was tested by experienced cardiac specialists in Germany.

The manufacture of your CVS blood pressure monitor is in accordance with the terms of the European standard for blood pressure measuring devices (see technical data) under the supervision of the Technical Monitoring Association Esses (RWTüV-Essen).

**Never attempt to repair the instrument yourself!**

Any unauthorized opening of the instrument invalidates all warranty claims!

## 8. Warranty

---

Your CVS blood pressure monitor is **guaranteed for a lifetime** against manufacturers' defects for the original purchaser only, from date of purchase. The warranty does not apply to damage caused by improper handling, accidents, professional use, not following the operating instructions or alterations made to the instrument by third parties.

There are no user serviceable parts inside. Batteries or damage from old batteries is not covered by the warranty. The warranty is only valid upon completion and delivery of the warranty card to address on the enclosed warranty card.

Please call toll free 1-866-464-6184 or email [CVSbpsupport@microlifeusa.com](mailto:CVSbpsupport@microlifeusa.com) for warranty claims.

**Please note: According to international standards, your monitor should be checked for calibration every 2 years.**

## 9. Certifications

---

### **Device standard:**

This device is manufactured to meet the European and United States standards for non-invasive blood pressure monitors  
EN1060-1 / 1995  
EN1060-3 / 1997

DIN 58130, NIBP – clinical investigation  
ANSI / AAMI SP10, NIBP – requirements

### **Electromagnetic compatibility:**

Device fulfills the stipulations of the European standard EN 60601-1-2

### **Clinical testing:**

Clinical performance tests were carried out in the US and Germany according to the DIN 58130/1997 procedure N6 (sequential) and AAMI standard (US).

The stipulations of the EU-Guidelines 93/42/EWG for Medical Products Class IIa have been fulfilled.

## 10. Technical specifications:

---

<b>Weight:</b>	465 g (with batteries)
<b>Size:</b>	131 (W) x 174 (L) x 73 (H) mm
<b>Storage temperature:</b>	-5 to +50°C
<b>Humidity:</b>	15 to 85% relative humidity maximum
<b>Operation temperature:</b>	10 to 40°C
<b>Display:</b>	LCD-Display (Liquid Crystal Display)
<b>Measuring method:</b>	oscillometric
<b>Pressure sensor:</b>	capacitive
<b>Measuring range:</b>	
<b>SYS/DIA:</b>	30 to 280 mmHg
<b>Pulse:</b>	40 to 200 per minute
<b>Cuff pressure display range:</b>	0-299 mmHg
<b>Memory:</b>	Stores last 30 measurements automatically
<b>Measuring resolution:</b>	1 mmHg
<b>Accuracy:</b>	Pressure within $\pm 3$ mmHg Pulse $\pm 5$ % of the reading
<b>Power source:</b>	a) 4 dry cells (batteries) UM-3, size AA, 1.5 V b) AC adaptor 6 V DC 600 mA (voltage 4.5 V DC to 6 V DC)
<b>Accessories:</b>	regular cuff for arm circumference 9.5"-13.25" large cuff for arm circumference 12"-16"

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©2004 CVS/pharmacy
**CVS Pharmacy, Inc.**